



ISQ and TSQ GC-MS

Mass Spectrometers

Spare Parts Guide

1R120621-0004 Revision A March 2022



© 2022 Thermo Fisher Scientific Inc. All rights reserved.

ISQ, TSQ, TRACE, and TriPlus are trademarks and Xcalibur is a registered trademark of Thermo Fisher Scientific in the United States.

The following are registered trademarks in the United States and other countries: Microsoft, Windows, Office, XP, and Excel are registered trademarks of Microsoft. Adobe, and Acrobat are registered trademarks of Adobe Systems Incorporated. SilTite is a registered trademark of SGE Analytical Science in the United States. Vespel is a registered trademark of E.I. du Pont de Nemours and Company TORX[®] is a registered trademark of Camcar LLC of Acument Global Technologies.

All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

Thermo Fisher Scientific Inc. makes no representations that this document is complete, accurate or errorfree and assumes no responsibility and will not be liable for any errors, omissions, damage or loss that might result from any use of this document, even if the information in the document is followed properly.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Revision A, March 2022.

General Lab Equipment, Not for Clinical, Patient or Diagnostic Use.

Contents

	Preface
	About Your System
	Related Documentationxii
	System Requirements
	Safety and Special Notices
	Special Notices
	Safety Symbols and Signal Words
	Hydrogen Safety Precautionsxv
	Using Hydrogen with the ISQ Series Mass Spectrometerxvi
	Using Hydrogen with TSQ Series Mass Spectrometer
	Hydrogen Connection Guidelinesxviii
	Purchasing Hydrogenxix
	Properly Storing Hydrogen
	Hydrogen Safety Codes, Standards and References
	Hazardous Substances Precautions
	Biological Hazard Warning Note
	Venting Toxic Gases
	Contacting Us
Chapter 1	Ordering Spare Parts
	Identifying A Part
	Calibration Gas Components
	Collision Gas Components
	Column Components
	Fan and Filter Components
	Ion Source Components
	Ion Source Cartridge Components
	Electron Multiplier Components
	Dual Filament Components
	Analyzer Components
	Analyzer Components
	Board Components
	Board Components
	Board Components
	Board Components

Pump Components	124
MS System Tools	
Upgrade Equipment	

Declaration

Manufacturer: Thermo Fisher Scientific

Thermo Fisher Scientific is the manufacturer of the instrument described in this manual and, as such, is responsible for the instrument safety, reliability and performance only if:

- installation,
- recalibration, and
- changes and repairs

have been carried out by authorized personnel and if:

- the local installation complies with local law regulations,
- the instrument is used according to the instructions provided, and
- if its operation is only entrusted to qualified trained personnel.

Thermo Fisher Scientific is not liable for any damages derived from the non-compliance with the aforementioned recommendations.

Regulatory Compliance

Thermo Fisher Scientific performs complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations. When the system is delivered to you, it meets all pertinent electromagnetic compatibility (EMC) and safety standards as described in the next section or sections by product name.

Changes that you make to your system may void compliance with one or more of these EMC and safety standards. Changes to your system include replacing a part or adding components, options, or peripherals not specifically authorized and qualified by Thermo Fisher Scientific. To ensure continued compliance with EMC and safety standards, replacement parts and additional components, options, and peripherals must be ordered from Thermo Fisher Scientific or one of its authorized representatives.

EMC and Safety Standards

- ITQ and Ion Trap Series standards: EMC EN 61326-1:2006. Safety IEC 61010-1:2001, IEC 61010-2-081:2001
- Direct Probe Controller (DPC) standards: EMC EN 61326-1:2006. Safety IEC 61010-1:2001, IEC 61010-2-081:2001
- ISQ, ISQ 7000, ISQ 7610 standards: EMC EN 61326-1:2013. Safety IEC 61010-1:2010, IEC 61010-2-010:2014, IEC 61010-2-081:2015.
- TSQ 8000, TSQ 8000 Evo, TSQ Duo, TSQ 9000, and TSQ 9610 standards: EMC EN 61326-1:2013. Safety IEC 61010-1:2010, IEC 61010-2-010:2014, IEC 61010-2-081:2015.

FCC Compliance Statement

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES. OPERATION IS SUBJECT TO THE FOLLOWING TWO CONDITIONS: (1) THIS DEVICE MAY NOT CAUSE HARMFUL INTERFERENCE, AND (2) THIS DEVICE MUST ACCEPT ANY INTERFERENCE RECEIVED, INCLUDING INTERFERENCE THAT MAY CAUSE UNDESIRED OPERATION.



CAUTION Read and understand the various precautionary notes, signs, and symbols contained inside this manual pertaining to the safe use and operation of this product before using the device.

Notice on Lifting and Handling of Thermo Scientific Instruments

For your safety, and in compliance with international regulations, the physical handling of this Thermo Fisher Scientific instrument *requires a team effort* to lift and/or move the instrument. This instrument is too heavy and/or bulky for one person alone to handle safely.

Notice on the Proper Use of Thermo Scientific Instruments

In compliance with international regulations: Use of this instrument in a manner not specified by Thermo Fisher Scientific could impair any protection provided by the instrument.

Notice on the Susceptibility to Electromagnetic Transmissions

Your instrument is designed to work in a controlled electromagnetic environment. Do not use radio frequency transmitters, such as mobile phones, in close proximity to the instrument.

EU REACH Statement

he European Commission promulgated legislation that covers the registration, evaluation, authorization and restriction of chemicals within the European Union community under (EC) No 1907/2006. This regulation is commonly known as REACH. Thermo Fisher Scientific is committed to meeting all compliance obligations under REACH. As per Article 33 of the Regulation, this product may include items which contain more than 0.1% by weight of some SVHC Candidate Substance. Some electronic parts and copper alloys can contain lead.

For manufacturing location, see the label on the instrument.

WEEE Directive

2012/19/EU

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with the following symbol:



Thermo Fisher Scientific is registered with B2B Compliance (B2Bcompliance.org.uk) in the UK and with the European Recycling Platform (ERP-recycling.org) in all other countries of the European Union and in Norway.

If this product is located in Europe and you want to participate in the Thermo Fisher Scientific Business-to-Business (B2B) Recycling Program, send an email request to weee.recycle@thermofisher.com with the following information:

- WEEE product class
- Name of the manufacturer or distributor (where you purchased the product)
- Number of product pieces, and the estimated total weight and volume
- Pick-up address and contact person (include contact information)
- Appropriate pick-up time
- Declaration of decontamination, stating that all hazardous fluids or material have been removed from the product

For additional information about the Restriction on Hazardous Substances (RoHS) Directive for the European Union, search for RoHS on the Thermo Fisher Scientific European language websites.

IMPORTANT This recycling program is not for biological hazard products or for products that have been medically contaminated. You must treat these types of products as biohazard waste and dispose of them in accordance with your local regulations.

Declarations of Conformity

		-Original-				
UK Declaratio	on of Conformity	UK CA				
SCIE	Thermo Fisher SCIENTIFIC					
2215 Grand	Fisher Scientific J Avenue Parkway Texas 78728 USA					
Declares, under sole responsibility, that	products					
Designation:	Mass Spectrometer					
Model:	Thermo Scientific TSQ Series TSQ9K Series, TSQ 9610 Series					
as originally delivered complies with the applicable UK Regulations:	essential requirements of the follo	wing				
Electrical Equipment (Safety) Regulations	2016					
Electromagnetic Compatibility Regulations	2016					
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (ROHS) Regulations	2012					
and complies with the following harmoni specifications:	zed standards and other technical					
BS EN 61010-1:2010+A1:2019	BS EN 61326-1:2021					
Signed for and on behalf of: Thermo Fisher Scientific:						
Brody Guckenberger (Director, Applied Research) Thermo Fisher Scientific						
13/3n	Austin, March 4, 2022					
Signature	Date					

-Original-

CE

EU Declaration of Conformity EU-Konformitätserklärung



Thermo Fisher Scientific 2215 Grand Avenue Parkway Austin, Texas 78728 USA

We hereby declare that the following products

Designation:	Mass Spectrometer

Thermo Scientific TSQ Series TSQ9K Series, TSQ 9610 Series

fulfill all the relevant requirements of the following directives:

Low Voltage Directive

2014/35/EU

Electromagnetic Compatibility Directive 2014/30/EU

RoHS Directive

Model:

2011/65/EU and (EU) 2015/863

The following relevant harmonized standards were used:

EN 61010-1:2020-03

EN 61326-1:2013-07

Person authorized to compile the technical file:

Brody Guckenberger (Director, Applied Research) Thermo Fisher Scientific

l

Signature

Austin, March 4, 2022 Date

Preface

This guide contains information about ordering spare parts for the Thermo ScientificTM ISQTM and TSQTM GC-MS mass spectrometers.

Contents

- About Your System
- Related Documentation
- System Requirements
- Safety and Special Notices
- Hydrogen Safety Precautions
- Hazardous Substances Precautions
- Contacting Us

About Your System

Gas chromatography/mass spectrometry (GC/MS) represents a combination of two powerful analytical techniques: GC, which acts as a separation technique, and MS, which acts as a detection technique. Complex mixtures of individual compounds can be injected into the GC, either manually or by an autosampler and then separated for presentation to the MS. The MS will generate a mass spectrum of the GC eluate and its components. The mass spectrum can then be used for qualitative identification as well as accurate and precise quantification of the individual compounds present in the sample.



WARNING Thermo Fisher Scientific systems operate safely and reliably under carefully controlled environmental conditions. If the equipment is used in manner not specified by the manufacturer, the protections provided by the equipment may be impaired. If you maintain a system outside the specifications listed in this guide, failures of many types, including personal injury or death, may occur. The repair of instrument failures caused by operation in a manner not specified by the manufacturer is specifically excluded from the Standard Warranty and service contract coverage.

Related Documentation

Thermo Fisher Scientific provides the following documents for the ISQ 7000 mass spectrometer:

- ISQ 7610 Preinstallation Requirements Guide, PN 1R120621-0001
- ISQ 7610 Hardware Manual, PN 1R120621-0002
- ISQ 7610 User Guide, PN 1R120621-0003
- ISQ and TSQ GC-MS Spare Parts Guide, PN 1R120621-0004
- Direct Probe System User Guide PN 1R120505-0006

Thermo Fisher Scientific provides the following documents for the TSQ 9000 mass spectrometer:

- TSQ 9610 Preinstallation Guide, PN 1R120622-0001
- TSQ 9610 User Guide, PN 1R120622-0002
- TSQ 9610 Hardware Manual, PN 1R120622-0003
- ISQ and TSQ GC-MS Spare Parts Guide, PN 1R120621-0004
- Direct Probe System User Guide, PN 1R120505-0006

Instrument Help is available from within the *AutoTune* and *Method Editor* software for each instrument.

System Requirements

The data system used with your ISQ 7610 or TSQ 9610 system must meet these minimum requirements:

System	Requirements
Hardware	 3.6 GHz dual-core processor enabled 16 GB RAM with system managed memory enabled DVD drive Resolution display 1280×1024 (SXGA) 20 GB available on drive C NTFS format
Software	 Microsoft[™] Windows[™] 10 Operating System (64-bit) English only or Windows 7 Professional Operating System (64-bit) Microsoft .NET Framework 4.0 or later Thermo Scientific[™] Xcalibur[™] and Foundation software^a Thermo Scientific[™] TraceFinder[™] software^b Thermo Scientific[™] Dionex[™] Chromeleon software^b

^a Check release notes for compatibility with ISQ Series or TSQ Series instrument control software.

^b Check release notes for compatibility with Thermo Foundation, Xcalibur, and SQ Series or TSQ Series instrument control software.

Safety and Special Notices

Make sure you follow the precautionary statements presented in this guide. The safety and other special notices appear in boxes.

Special Notices

Special notices include the following:

IMPORTANT Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

Note Highlights information of general interest.

Tip Highlights helpful information that can make a task easier.

Safety Symbols and Signal Words

All safety symbols are followed by **WARNING** or **CAUTION**, which indicates the degree of risk for personal injury, instrument damage, or both. Cautions and warnings are following by a descriptor. A **WARNING** is intended to prevent improper actions that *could* cause personal injury. A **CAUTION** is intended to prevent improper actions that *might* cause personal injury or instrument damage. You can find the following safety symbols on your instrument or in this guide.

Symbol	Descriptor
	BIOHAZARD: Indicates that a biohazard <i>will, could</i> , or <i>might</i> occur.
	BURN HAZARD: Alerts you to the presence of a hot surface that <i>could</i> or <i>might</i> cause burn injuries.
	ELECTRICAL SHOCK HAZARD: Indicates that an electrical shock <i>could</i> or <i>might</i> occur.

Symbol	Descriptor
	FIRE HAZARD: Indicates a risk of fire or flammability <i>could</i> or <i>might</i> occur.
RAMMABLE GAS 2	FLAMMABLE GAS HAZARD: Alerts you to gases that are compressed, liquefied or dissolved under pressure and can ignite on contact with an ignition source. This symbol indicates this risk <i>could</i> or <i>might</i> cause physical injury.
	GLOVES REQUIRED: Indicates that you must wear gloves when performing a task or physical injury <i>could</i> or <i>might</i> occur.
	HAND AND CHEMICAL HAZARD: Indicates that chemical damage or physical injury <i>could</i> or <i>might</i> occur.
	INSTRUMENT DAMAGE: Indicates that damage to the instrument or component <i>might</i> occur. This damage might not be covered under the standard warranty.
	LIFTING HAZARD: Indicates that a physical injury <i>could</i> or <i>might</i> occur if two or more people do not lift an object.
	MATERIAL AND EYE HAZARD: Indicates that eye damage <i>could</i> or <i>might</i> occur.
	RADIOACTIVE HAZARD: Indicates that exposure to radioactive material <i>could</i> or <i>might</i> occur.
2	READ MANUAL: Alerts you to carefully read your instrument's documentation to ensure your safety and the instrument's operational ability. Failing to carefully read the documentation <i>could</i> or <i>might</i> put you at risk for a physical injury.

Symbol Descriptor



For the prevention of personal injury, this general warning symbol precedes the **WARNING** safety alert word and meets the ISO 3864-2 standard. In the vocabulary of ANSI Z535 signs, this symbol indicates a possible personal injury hazard exists if the instrument is improperly used or if unsafe actions occur. This symbol and another appropriate safety symbol alerts you to an imminent or potential hazard that *could cause personal injury*.

TOXIC SUBSTANCES HAZARD: Indicates that exposure to a toxic substance could occur and that exposure *could* or *might* cause personal injury or death.

Hydrogen Safety Precautions

Hydrogen is a colorless, odorless, highly flammable gas with the molecular formula H_2 . Hydrogen gas presents a hazard as it is combustible over a wide range of concentrations: at ambient temperature and pressure, this ranges from about 4% to 74.2% by volume.

Hydrogen has a flash point of - 423 °F (- 253 °C) and an auto-ignition temperature of 1,040 °F (560 °C). It has a very low ignition energy and the highest burning velocity of any gas. If hydrogen is allowed to expand rapidly from high pressure, it can self-ignite. Hydrogen burns with a flame that can be invisible in bright light.



WARNING FIRE HAZARD: The use of hydrogen as a carrier gas is dangerous. Hydrogen is potentially explosive and must be used with extreme care. Any use of hydrogen gas must be reviewed by appropriate health and safety staff and all installations of hydrogen systems must be performed to applicable codes and standards. Thermo Fisher Scientific assumes no liability for the improper use of hydrogen as a carrier gas.

Before you begin using hydrogen, you should conduct a risk assessment based on the quantity of hydrogen to be used and the conditions of your laboratory. You should ask yourself:

"What hydrogen hazards associated with this project are most likely to occur?"

"What hydrogen hazards associated with this project have the potential to result in the worst consequences?"

- Try to reduce or eliminate the higher risks by using the proper ventilation to remove hydrogen gas before an ignitable concentration can accumulate. You should also consider purging the hydrogen to further reduce hazards and ensure anyone who will be working with hydrogen has basic hydrogen safety training.
- As with laboratory safety in general, be sure to wear safety glasses, laboratory coats, gloves, etc. Typically there are no specific requirements for gaseous hydrogen, other than eye protection when working with a compressed gas. If working with liquid (cryogenic) hydrogen, insulated gloves and protective shoes should be worn in addition to eye protection.

- You should post "No Smoking" and "No Open Flames" signs to identify hydrogen sources and cylinders. Maintain, inspect and leak-test all hydrogen sources regularly.
- All hydrogen shutoff valves should be clearly marked and permanent hydrogen piping should be labeled as such at the supply or discharge point and at regular intervals along its length. Where hydrogen gas piping passes through a wall, the piping should be labeled on both sides of the wall.
- There should also be contingency plans in place should an incident occur.
- The site emergency response team, as well as the local fire department, should know the location of all hydrogen storage tanks.

Using Hydrogen with the ISQ Series Mass Spectrometer

To use hydrogen with the ISQ Series instrument, you must always shut off the GC carrier gas before venting or turning off the ISQ Series instrument. There are three hydrogen safety screws on the ISQ Series instrument that **must** be in place. These are attached to your instrument at the factory.




Figure 1. Hydrogen Safety Screws on the ISO Series Instrument

Before powering on the ISQ Series system, ensure that:

- All the covers and panels of the ISQ Series system are firmly attached.
- The vent valve is tightly closed if you vented the system.
- All fittings, ferrules, and o-rings are sealed.

Using Hydrogen with TSQ Series Mass Spectrometer

To use hydrogen with the TSQ Series instrument, you must always shut off the GC carrier gas before venting or turning off the TSQ Series instrument. There are three hydrogen safety screws on the TSQ Series instrument that **must** be in place. These are attached to your instrument at the factory.




Figure 2. Hydrogen Safety Screws on the TSQ Series Instrument

Before powering on the TSQ Series system, ensure that:

- All the covers and panels of the TSQ Series system are firmly attached.
- The vent valve is tightly closed if you vented the system.
- All fittings, ferrules, and o-rings are sealed.

Hydrogen Connection Guidelines

Use the following guidelines to safely connect hydrogen to your system:

• **Piping**—Hydrogen must be delivered to equipment using appropriate piping and be done in such a way as to pose essentially no hazard to end-users. Piping systems for the delivery of hydrogen should be designed and installed by a person qualified by specific training and experience with hydrogen piping systems.

Stainless steel is usually recommended because it is a safe, cost-effective material. Piping of *black iron* or copper must not be used, as the pipe can become brittle with age. Elastomeric/plastic tubing of various plastics and polymers should not be used, unless the tubing is approved for use with hydrogen. If elastomeric/plastic tubing is used for hydrogen gas delivery, the tubing should be tested for hydrogen permeability to minimize leakage.

The hydrogen piping system must be flexible enough to endure routine thermal expansion and contraction. The system should also include considerations for the most severe condition of temperature and pressure expected during service. Piping and supports must be able to withstand static loading introduced by such things as ice and snow; and dynamic loading from high wind and earthquake.

Caution should be used if burying hydrogen piping. Proper controls should be used to protect against damage and corrosion, and also to prevent Hydrogen from entering a building if there is any leakage.

• Fittings—All fittings must be of the proper type approved or designed for use with hydrogen gas. Use as few fittings as possible to minimize the potential for leaks. After installation, ensure that leak testing is carried out prior to system use, and on a regular basis.

There must be no PTFE tape or other things like *plumber's putty* used to enhance a seal, as this actually is a detriment to a good seal. Ideally the best installation would use stainless steel tubing with appropriate gas-tight fittings.

Welding is usually preferred for joints in hydrogen piping systems since welding provides a better connection and reduces the potential for leaks compared to mechanical fittings. Soft solder joints are not permitted for hydrogen systems (due to the low melting point of soft solder and its potential for brittle failure at cryogenic temperatures). Brazed joints are permitted, but such joints should be protected against the possibility of external fire.

Tubing connections should be clamped to barbed or press-fit type connections. Hose clamps or *jubilee clamps* must not be used.

• Valves—All valves must be suitable for hydrogen service and for the specific operating conditions. Valves, including regulators, must not be used for hydrogen, unless they are designed and identified for such a use. Ball valves are often chosen because of their superior leak tightness through the valve seat. Pneumatic operators are usually chosen for remotely operated valves so that potential ignition sources (electricity) are remote from the valve.

Manual shutoff valves should be provided near each point of use, within immediate reach. If a hydrogen cylinder or hydrogen generation system is located within immediate reach, a separate point-of-use shutoff valve is usually not necessary.

Line regulators that have their source away from the point of use should have a manual shutoff valve near the point of use.

An emergency gas shutoff device in an accessible location outside the use area should be provided in addition to the manual point-of-use valve in each educational and instructional laboratory space that has a piped gas supply system.

If necessary, the piping system should have uninterruptible pressure relief. The pressure relief system should be designed to provide a discharge rate sufficient to avoid further pressure increase and should vent to a safe location outside or to a ventilation system exhaust.

Purchasing Hydrogen

Use the following guidelines when purchasing hydrogen:

• Hydrogen Generator—Because it minimizes the amount of hydrogen present and reduces the degree of hazard, a hydrogen generator (also called an electrolyzer) is the safest way to purchase hydrogen in the quantity used in GC/MS.

However, to minimize the degree of hazard, the hydrogen generator must only be operated in a non-explosive environment because hydrogen buildup can be ignitable. This means that your ventilation system for the room or lab hood must maintain an air exchange rate that is at least two orders of magnitude greater than the maximum hydrogen production rate of the hydrogen generator. Be sure to follow the manufacturers' directions about proper use and maintenance of the regulator.

To prevent the possibility of releasing hydrogen, the hydrogen generator should be set to shut down if:

- There is a loss of flow to the ventilation system
- A hydrogen detector alarms at 25% of the lower flammable limit of hydrogen in air.

The oxygen exhausted by the electrolyzer should be vented to the outside as well.

• Hydrogen Cylinder—Hydrogen can be delivered in standard laboratory gas bottles or cylinders. These cylinders have a limited amount of hydrogen in them and are a safe way to transport and store hydrogen. However, compressed hydrogen gas cylinders, like all compressed gas cylinders, must be secured in an upright position, ideally with a non-combustible chain or cable. If the cylinder falls over, the valve can be knocked off and the pressurized cylinder can take off like a rocket, which leads to the release of hydrogen and possibly an explosion, severe injury, or death. Never crack a hydrogen cylinder valve to remove dust or dirt from fittings prior to attaching a regulator, as there is a risk of self-ignition.

Properly Storing Hydrogen

Storing and handling compressed hydrogen gas and cryogenic liquid hydrogen present potential health and safety hazards. Using proper storage and handling techniques is essential to maintaining a safe work environment.

Use the following guidelines when storing hydrogen:

- Store spare hydrogen gas cylinders outside and away from doors, windows, building air intake vents, structures, and vehicle routes. This precaution applies when the hydrogen is or is not in use. Indoor storage of spare hydrogen cylinders has special requirements, which is beyond the scope of this document. Documentation for each vessel should include a description of the vessel, a list of available drawings or other documents, the most recent inspection results, and the responsible person's name.
- Prevent spare cylinders from toppling by wrapping them with chains. The chains should also be protected against corrosion and excessive heat.
- Separate spare hydrogen cylinders from oxidizing gases (such as oxygen) with a 5 ft (1.5 m) tall fire barrier with a half-hour fire rating or place the cylinders at least 20 ft (6 m) apart.
- When moving hydrogen cylinders:
 - Remove the regulator and replace the cylinder valve cap before moving.
 - Move cylinders on cylinder carts or with other appropriate transport devices.
 - Never roll or drop a cylinder and never lift a cylinder by its protective cap.
- Bulk hydrogen systems include either gaseous or liquid hydrogen in fixed installations; in some gas systems a semi-permanent trailer (tube trailer) can be used. Storage vessels for compressed hydrogen gas or liquid hydrogen should be designed, constructed, tested, and maintained in accordance with applicable codes and standards. Bulk hydrogen systems represent a level of complexity again which is beyond the scope of this document; however some general guidelines are provided.
- The bulk hydrogen storage system should not be located beneath electric power lines, close to other flammable gases/liquids, or close to public areas. It should be readily accessible to authorized personnel and delivery equipment, but protected from physical damage or tampering.
- As liquid hydrogen systems also have a cryogenic hazard, additional safety considerations for the use of cryogenic liquids may be necessary.

Hydrogen Safety Codes, Standards and References

The following list of safety codes, standards and references is in no way an exhaustive list. In fact, there may be federal, state or local codes that apply to your specific location. Check with all appropriate agencies with jurisdiction before installing or using a hydrogen system.

- Air Products Safetygram #4 Gaseous Hydrogen
- ANSI/AIAA standard for hydrogen safety guidelines is AIAA G-095-2004, Guide to Safety of Hydrogen and Hydrogen Systems
- ASME B31.1, Power Piping Code
- ASME B31.3, Process Piping Code
- ASME B31.8, Gas Transmission and Distribution Systems
- BCGA Code Of Practice CP4 Industrial Gas Cylinder Manifolds and Gas Distribution Pipework
- BCGA Code Of Practice CP33 The Bulk Storage of Gaseous Hydrogen at Users' Premises
- CGA G-5, Hydrogen
- CGA G-5.4, Standard for Hydrogen Piping Systems at Consumer Locations
- CGA G-5.5, Hydrogen Vent Systems
- CGA G-5.6, Hydrogen Pipeline Systems
- CGA G-5.8, High Pressure Hydrogen Piping Systems at Consumer Locations.
- FM Global Property Loss Prevention Data Sheets 7-50: Compressed Gases in Cylinders
- FM Global Property Loss Prevention Data Sheets 7-91: Hydrogen
- IGC Doc 121/04/E, Hydrogen Transportation Pipelines System Design Features
- NASA
- NSS 1740.16 Safety Standard For Hydrogen And Hydrogen Systems Guidelines for Hydrogen System Design, Materials Selection, Operations, Storage, and Transportation
- NFPA 52, Vehicular Fuel Systems Code
- NFPA 55, Standard for the Storage, Use, and Handling of Compressed Gases and Cryogenic Fluids in Portable and Stationary Containers, Cylinders, and Tanks, 2005 Edition
- NFPA 68, Standard on Explosion Protection by Deflagration Venting
- NFPA 70, National Electrical Code

- NFPA 497, Recommended Practice for the Classification of Flammable Liquids, Gases, or Vapors and of Hazardous (Classified) Locations for Electrical Installations in Chemical Process Areas
- NFPA 13, Standard for the Installation of Sprinkler Systems
- NFPA 45, Standard on Fire Protection for Laboratories Using Chemicals
- NFPA 55, Standard for the Storage, Use, and Handling of Compressed Gases and Cryogenic Fluids in Portable and Stationary Containers, Cylinders, and Tanks
- NFPA 68, 2007 Standard on Explosion Protection by Deflagration Venting
- NFPA 69, Standard on Explosion Prevention Systems
- NFPA 91, Standard for Exhaust Systems for Air Conveying of Vapors
- NFPA 255, Standard Method of Test of Surface Burning Characteristics of Building Materials
- OSHA 29CFR1910.103 1910.103 Hydrogen

Hazardous Substances Precautions



WARNING Before using hazardous substances (toxic, harmful, and so on), please read the hazard indications and information reported in the applicable Material Safety Data Sheet (MSDS). Use personal protective equipment according to the safety requirements.

Biological Hazard Warning Note

In laboratories where samples with potential biological hazards are handled, the user must label any equipment or parts which might become contaminated with biohazardous material.



The appropriate warning labels are included with the shipment of the instrument. It is the user's responsibility to label the relevant parts of the equipment.

When working with biohazardous materials, you are responsible for fulfilling the following mandatory requirements:

- Providing instructions on how to safely handle biohazardous material.
- Training operators to be aware of potential hazards.

- Providing personal protective equipment.
- Providing instructions for what to do if operators are exposed to aerosols or vapors during normal operation (within the intended use of the equipment) or in case of single fault situations such as a broken vial. The protective measures must consider potential contact with the skin, mouth, nose (respiratory organs), and eyes.
- Providing instructions for decontamination and safe disposal of relevant parts.



WARNING The user or operator is responsible for the safe handling of hazardous chemicals or biological compounds including (but not limited to) bacterial or viral samples and the associated waste, according to international and local regulations.

Venting Toxic Gases

When analyzing toxic compounds, be aware that during the normal operation of the GC some of the sample might be vented outside the instrument through the split and purge flow vents; therefore, be sure to vent the exhaust gases to a fume hood. Consult local environmental and safety regulations for instructions in exhausting fumes from your system.

Contacting Us

There are several ways to contact Thermo Fisher Scientific for the information you need.

To find out more about our products

Go to www.thermofisher.com for information about our products.

✤ To get local contact information for sales or service

Go to gcgcms.freshdesk.com/support/home

- To suggest changes to documentation or to Help
 - Fill out a reader survey online at www.surveymonkey.com/s/PQM6P62.
 - Send an e-mail message to the Technical Publications Editor at techpubs-austin@thermofisher.com.

Ordering Spare Parts

This chapter contains illustrations and part numbers for all of the replaceable components in ISQ and TSQ GC-MS instruments. Refer to the hardware manual for your instrument for information about installing these components in your instrument. Any component throughout the documentation that has a part number can be ordered from us. Components without part numbers are not available.

Contents

- Identifying A Part
- Calibration Gas Components
- Collision Gas Components
- Column Components
- Fan and Filter Components
- Electron Multiplier Components
- Dual Filament Components
- Analyzer Components
- Board Components
- Power Supply Components
- Vacuum Interlock Components
- Manifold Components
- Cover Components
- Pump Components
- MS System Tools
- Upgrade Equipment

Identifying A Part

To identify a part, you need to know where it is located in the instrument or the part's relationship to a particular functionality of the instrument. Use the categories in the table below to find a location or functionality that relates to the component you need. For example, if you know the part is related to the calibration gas, look in the Calibration Gas Components category. Then refer to the illustrations in that category to visually identify that particular part.

IMPORTANT Only components with a part number are available for purchase. The part you are looking for may be shown in an illustration, but if it does not list a part number, it is not available.

Table 1. Parts That Can Be Reordered

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Single Level Calibration Gas Module	1R119246-0004	All	
Dual Level Calibration Gas Module	1R119246-0003	All	
Calibrant Reservoir Kit	1R120433-0001	All	
Graphite/Vespel Ferrule for Calibration Gas Modules	1R76458-0030	ISQ 7000 and ISQ 7610; TSQ 9000 and TSQ 9610	
Transfer Line	1R120402-0001	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ Duo, TSQ 8000 Evo	
O-ring for Transfer Line	1R3814-223	All	
Screws for Transfer Line	1R76913-0410	All	
Gas Mixing Chamber Kit	1R120574-KIT	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ Duo, TSQ 8000 Evo	
Gas Mixing Chamber (1)	1R1205740054	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ Duo, TSQ 8000 Evo	Calibration Gas Components
Source Gas Tube (2 in Kit)	1R120574-0055 (PN orders 1)	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ Duo, TSQ 8000 Evo	
M3 × 8 mm Screw (2 in Kit)	1R76913-0308	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ Duo, TSQ 8000 Evo	
Calibration Gas Line	1R120438-0010	All	
Vacuum Hose (Order 10 ft)	1R76505-0625	All	
FC-43 Calibration Compound	1R50010-30059	All	
1 fg/μL OFN IDL Kit	1R76310-0101K	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	-
5 fg/µL OFN in Iso-octane	1R76310-0105	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	-
ISQ 7000, ISQ 7610, TSQ 9000 and TSQ 9610 Transfer Line Assembly	1R120191-0020	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
Calibration Gas Tube	1R120191-3000	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Collision Gas Module	1R119246-0005	TSQ All	
Tube Fitting, Collision Gas Module to Manifold	1R76256-0045	TSQ All	Collision Gas
Collision Gas Tube	1R120557-0030	TSQ All	
Tee, 1/8 in.	1R76256-0126	TSQ All	
Graphite Vespel [™] Ferrule for 0.18 mm Columns (pkg of 10)	1R76458-2016	All	
Graphite Vespel Ferrule for 0.25 mm Column (pkg of 10)	29033496	All	
Graphite Vespel Ferrule for 0.32 mm Columns (pkg of 10)	1R76458-2019	All	
Graphite Vespel Ferrule for 0.53 mm Columns (pkg of 10)	1R76458-2020	All	
Graphite Vespel Ferrule for Calibration Gas Module Connection	1R176458-0030	All	
2-hole Graphite Vespel Ferrule for <0.32mm Column (pkg of 10)	1R76458-2018	All	
No-Hole Graphite Vespel Ferrule (pkg of 10)	1R76458-2009	All	Column Components
Graphite Vespel Ferrule for 0.1–0.25 mm Columns to be used with Spring-loaded Transfer Line Nut (pkg of 10)	290VA191	All	
Spring-loaded Transfer Line Nut	1R120434-0010	All	
Nickel-Coated Nut for Graphite Vespel Ferrule (pkg of 5)	1R76256-0060	All	
SilTite Ferrule for 0.10-0.25mm ID Column (pkg of 10)	1R76458-2000	All	
SilTite Ferrule for 0.32mm Columns (pkg10)	1R76458-2024	All	
SilTite Ferrule for 0.53mm Columns (pkg10)	1R76458-2026	All	
Nut for SilTite Ferrules (pkg of 5)	1R76458-2001	All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Intake Filter/RF Shield (1)	1R76475-5000	All	
Triple Filter with Base	A0950-01600	All	
Triple Filter Replacement Cartridge	A0950-R1600	All	
Oil Mist Filter	1R76505-0036	All	
Chassis Cooling Fan	1R120443-0001	All	
Nylon Rivet for Cooling Fan	1R3326-5000	All	Fan and Filter
Air Deflector	1R120552-0010	TSQ All	Components
Front Wire Mesh Filter	1R120554-0010	TSQ All	
Left Rear Wire Mesh Filter	1R120554-0020	TSQ All	
Fan Plenum	1R120553-0020	TSQ All	
Fan Filter (Plastic)	1R76475-5005	TSQ All	
Turbo Cooling Fan	1R120443-0003	TSQ All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Ion Source Block	1R120404-1311	ISQ, ISQ LT, ISQ QD, TSQ 8000, TSQ Duo, TSQ 8000 Evo	
ISQ 7000, ISQ 7610, TSQ 9000 Ion Source Block	1R120191-2000	ISQ 7000, ISQ 7610, TSQ 9000	
TSQ 9610 Ion Source Block	1R120189-2024	TSQ 9610 with NVAEI	
NVAEI Lock Plate Assembly	1R120189-6015	TSQ 9610 with NVAEI	
NVAEI Lock Plate Washer	1R120189-2060	TSQ 9610 with NVAEI	
NVAEI Heat Block Alignment Pin	1R120189-2022	TSQ 9610 with NVAEI	
NVAEI Source Exchange Tool	1R120189-0003	TSQ 9610 with NVAEI	
NVAEI Source Sleeve	1R120189-2000	TSQ 9610 with NVAEI	
NVAEI Ion Volume	1R120189-2001	TSQ 9610 with NVAEI	
NVAEI Source Insulator Filament	1R120189-2036	TSQ 9610 with NVAEI	
NVAEI Ion Source	1R120189-6030	TSQ 9610 with NVAEI	
NVAEI Dual Filament	1R120189-6000	TSQ 9610 with NVAEI	
Repeller Plate	1R120404-1312	All	
Source-Repeller Thumbscrew	1R120564-0004	ISQ, ISQ LT, ISQ QD, TSQ 8000, TSQ Duo, TSQ 8000 Evo	-
ISQ 7000, ISQ 7610, TSQ 9000, and TSQ 9610 Source-Repeller Thumbscrew	1R120602-0206	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
Thumbscrew-Repeller Spacer	1RTSLI053A	All	
Lens/Source Heater	1R120404-1330	All	Ion Source Components
Lens Heater Block	1R120404-1321	All	Components
Lock Washer for Lens Heater Block	1R76902-0004	All except TSQ 9610 with NVAEI	
Screw for Lens Heater Block	1R76913-0306	All	
Lens Heater Block Grounding Strap	1R120404-2010	All	1
Screw for Grounding Strap	1R76913-0305	All	1
Lens Plate and Springs	1R120404-1750	All	1
Copper Washer for Lens Plate	1R76483-5003	All	1
Screw for Lens Plate	1R76913-0306	All	1
Source to Ion Guide Spacer	1R120404-1340	All]

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
• Lens 1	1R120574-1039	All	
• Lens 2	1R120404-1140	All	
Electron Multiplier Plate	1R120405-1000	All	
Electron Multiplier	1R76022-14633	All except ISQ 7610 and TSQ 9610	-
Electron Multiplier, ISQ 7610 and TSQ 9610	1R76022-14634	ISQ 7610 and TSQ 9610	-
Dynode Feedthrough	1R120610-0040	All	
Detector Assembly (Dynode, Electron Multiplier, and Lead Set)	1R120379-0001	All	Electron Multiplier
O-ring for Large Feedthrough	1R3814-113	All	Components
Screw for Large Feedthrough	1R76913-0410	All	
Anode Feedthrough Assembly Kit, includes:	1R120480-0012	All	
Anode Feedthrough Assembly	1R120610-0050	All	
O-ring for Anode Feedthrough	1R3814-110	All	
Standoff for Anode Feedthrough	1R77005-3010	All	-
Screw for Anode Feedthrough (pkg of 4)	1R76913-0306	All	
Dual Filament	1R120404-1940	All	
Filament Retaining Spring Kit	1R120404-1405	All	Dual Filament Components
Screw for Filament Retaining Spring	1R76913-0306	All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
ISQ EI Analyzer Tray Assembled	1R120404-0003	ISQ, ISQ LT, ISQ QD	_
ISQ 7000 and ISQ 7610 Analyzer Tray Assembled	1R120191-0011T	ISQ 7000, ISQ 7610	
TSQ Analyzer Tray Assembled	1R120542-2535	TSQ 8000, TSQ 8000 Evo, TSQ Duo	-
TSQ 9000 Analyzer Tray Assembled	1R120192-0011T	TSQ 9000	
TSQ 9610 with NVAEI Analyzer Tray, Assembled	1R120189-0011T	TSQ 9610 with NVAEI	-
ISQ Analyzer Tray	1R120404-2001	ISQ All	
Ion Guide	1R120404-3100	All	
Ion Guide Clamp	1R120404-3214	All	
Screws for Ion Guide Clamp	1R76913-0316	All	Analyzer Components
Q1 Assembly with Wires and Tested	1R120542-2560	All	
Q1 Exit Endcap	1R120404-3230	ISQ All	
Q1 Entrance Endcap	1R120404-3211	All	
Q1 Quad Entrance Lens	1R120404-3212	All	
Quad 1 Wire Set	1R120480-0005	All	
Quad 3 Wire Set	1R120480-0006	TSQ All	
Q3 Quad Assembly	1R120542-2570	TSQ All	
Q1 Exit Lens Wire	1R120542-2515	TSQ All	1
Q3 Entrance Lens Wire	1R120542-2514	TSQ All	1
Q3 Exit Lens Wire	1R120542-2513	TSQ All	1
Q3 Exit Lens Assembly	1R120542-2536	TSQ All	1
Wire Routing Clip	1R120542-2511	TSQ All	1

Thermo Scientific Component Used on Instrument(s) Category Part Number ISQ All ISQ Distribution Board 1R120369-0001 TSQ Distribution Board 1R120485-0810 TSQ All Screw for Distribution Board 1R76913-0310 All Q1/Front RF Board 1R120485-A060 All Q3/Rear RF Board TSQ All 1R120485-A065 Q1/Front Rod Driver Board 1R120485-0415 All Q3/Rear Rod Driver Board (75V) 1R120485-0420 TSQ All Standoff for Rod Driver Board 1R76350-5003 All Screw for Rod Driver Board 1R76913-0310 All Lens Driver Board-ISQ 1R120354-A110 ISQ All Lens Driver Board-TSQ 1R120485-A110 TSQ 8000, TSQ Duo, TSQ 8000 Evo TSQ 9000 and TSQ 9610 Lens Driver 1R120485-A114T TSQ 9000, TSQ 9610 Board Controller Interface Board 1R120485-0020 ISQ All **Board Components** Controller Interface Board Support Bracket ISQ All 1R120373-0002 Source Interface Board, ISQ Systems 1R120354-0210 ISQ All Source Interface Board, TSQ Systems 1R120485-0210 TSQ All Heat Shield for AEI Source Interface Board TSQ 9610 with NVAEI 1R120189-2111 Electrometer Board 1R120354-0500 All Screw (6 mm) for Electrometer Board 1R76913-0306 All Screw (16 mm) for Electrometer Board 1R76913-0316 All Heat Shield for Source Interface Board All 1R120404-2110 Screw for Heat Shield 1R76913-0410 All 15-pin Male/Female RS-232 Cable All 1R76396-0500 Electrometer Shield All 1R120368-0010 Screw for Electrometer Shield All 1R76913-0306 O-ring for 4-Pin Feedthrough All 1R3814-127 All Screw for 4-Pin Feedthrough 1R76913-0410 20-pin Feedthrough 1R120610-0020 All

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
O-ring for 20-Pin Feedthrough	1R3814-123	All	
Screw for 20-Pin Feedthrough	1R76913-0410	All	
4-pin Feedthrough	1R120610-0030	All	-
Q1 Rod Driver Fan, 5 VDC	1R120443-0011	All	
Q1 RF Fan, 60 mm	1R120443-0002	All	-
Q3 RF Fan	1R120443-0005	All	-
PC Communication Board	1R120595-0010	All	-
ISQ Power Supply System	1R120380-0001	ISQ All	
TSQ Power Supply System	1R120544-0002	TSQ All	Power Supply
Dynode and Multiplier Power Supply and Cables (EI only)	1R120361-0003	All	Components
Inner Ball Seal	1R120406-1008	All	
Valve Axle	1R120406-1003	All	-
Axle Seal and O-ring	1R120406-2204	All	
Axle Bearing	1R120406-1005	All	-
Sealing Ball	1R120406-1004	All	-
Outer Ball Seal	1R120406-1002	All	
Handle	1R120406-1200	All except TSQ 9610 with AEI VPI	
Handle, VPI AEI Systems	1R120189-6020	TSQ 9610 with NVAEI	-
Solenoid Valve, 24 V	1R76461-5002	All	Vacuum Interlock
Barb Fitting, "L"	1R76256-1332	All	Components
Barb Fitting, Straight	1R76256-0040	All	
Valve Spring (pkg of 5)	1R76485-1001	All	
Vacuum Interlock Knob	1R120406-3000	All	
Vacuum Interlock Assembly	1R120403-0001	All	
Tubing (cm)	1R76433-0107	All	
Vacuum Interlock Microswitch Assembly	1R120406-1030	All	
Interior O-ring	1R3814-127	All	1
Exterior O-ring	1R3815-320	All	-
Clip (Pkg of 5)	1R76483-2102	All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Manifold Door Hinge	1R120565-0001	All	
Front Manifold Plate Assembly, with VPI	1R120403-0101	ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo	
Front Manifold Plate Assembly, with VPI	1R120191-0012	ISQ 7000 VPI, TSQ 9000 VPI	
Front Manifold Plate Assembly, with VPI, AEI	1R120189-0012	TSQ 9610 with NVAEI	-
Front Plate VPI AEI Adapter	1R120189-2056	TSQ 9610 with NVAEI	
Door Alignment Socket for VPI AEI	1R120189-2059	TSQ 9610 with NVAEI	
AEI Magnet Rotator	1R120189-6201	TSQ 9610 with NVAEI	
AEI Lock Plate Replacement	1R120189-6202	TSQ 9610 with NVAEI	
Front Manifold Plate Assembly, No VPI	1R120403-0002	ISQ QD, TSQ Duo	
Front Manifold Plate Assembly, No VPI	1R120191-0013	ISQ 7000 No VPI, TSQ 9000 No VPI	Manifold Components
O-ring for Front, Top, Detector, and Back Plate	1R3815-360	All	
Screw for Front Plate (M4 × 20 mm)	1R176913-0420	All	
Vent Valve Knob	1R120403-0104	All	
Source Insertion Guide	1R120564-0001	All	
O-ring for Vent Valve Knob	1R3814-110	All	
Back Tray Alignment Pin	1R120403-0103	All	
Front Door Alignment Pin	1R120564-0003	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo	
Source Magnet Retainer	1R120602-0205	ISQ 7000, TSQ 9000	-
Back Manifold Plate	1R120403-1001	All	
Magnet Yoke	1R120564-0002	All	
Source Magnet	1R70001-98195	ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo	

Right Side Panel-ISQ oldIR120413-0001ISQ, ISQ QD, ISQ ITRight Side Panel-ISQIR120191-1020ISQ 7000Left Side Panel-ISQIR120411-0001ISQ AllTop Cover Panel-ISQIR120412-0001ISQ AllTop Manifold Cover (glass)IR120191-0110ISQ 7610 with VP1ISQ 7610 Front Cover (with VP1)IR120191-0210ISQ 7610 no VP1ISQ 7610 Front Cover (with VP1)IR120191-0210ISQ 7000 no VP1ISQ 7000 Front Cover (with VP1)IR120191-0200ISQ 7000 with VP1ISQ 7000 Front Cover (with VP1)IR120191-0200ISQ 7000 with VP1ILarge Top Glass Cover-TSQIR12047-0200ISQ AllSught Side Panel-TSQIR120541-0000TSQ AllSught Side Panel-TSQIR120541-0001TSQ AllSught Side Panel-TSQIR120541-0003TSQ AllSught Side Panel-TSQIR120541-0001TSQ AllSught Side Panel-TSQIR120541-0003TSQ AllSug 610 Front Cover (with VP1)IR12019-0120TSQ Soft With VP1TSQ 9610 Front Cover (with VP1)IR12019-0120TSQ 9610 with VP1TSQ 9610 Front Cover (with AEI VP1)IR12019-0120TSQ 9600 Fon VP1TSQ 8000 Evo Front Cover (with AEI VP1)IR120548-0300TSQ 8000 Evo no VP1TSQ 9610 Front Cover (with VP1)IR120548-0300TSQ 8000 Evo no VP1TSQ 9610 Front Cover (with VP1)IR12019-0100TSQ 9000 Fon VP1TSQ 9000 Front Cover (with VP1)IR12019-0100TSQ 9000 Fon VP1TSQ 9000 Front Cover (with VP1)IR12019-0100TSQ 9000 with VP1	Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Left Side Panel-ISQIR120411-0001ISQ AllTop Cover Panel-ISQIR120412-0001ISQ AllTop Manifold Cover (glass)IR120401-3000ISQ AllISQ 7610 Front Cover (with VPI)IR120191-0110ISQ 7610 with VPIISQ 7610 Front Cover (no VPI)IR120191-0210ISQ 7610 no VPIISQ 7610 Front Cover (with VPI)IR120407-0100ISQ DISQ Front Cover (with VPI)IR120191-0210ISQ 7000 with VPIISQ 7000 Front Cover (with VPI)IR120191-0200ISQ 7000 no VPIISQ 7000 Front Cover (with VPI)IR120191-0200ISQ 7000 no VPIIArge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0200TSQ 9000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 no VPITSQ 9000 Front Cover	Right Side Panel-ISQ old	1R120413-0001	ISQ, ISQ QD, ISQ LT	
Top Cover Panel-ISQIR120412-0001ISQ AllTop Manifold Cover (glass)IR120401-3000ISQ AllISQ 7610 Front Cover (with VPI)IR120191-0110ISQ 7610 with VPIISQ 7610 Front Cover (no VPI)IR120191-0210ISQ 7610 no VPIISQ Front Cover (with VPI)IR120407-0100ISQ DISQ 7000 Front Cover (with VPI)IR120191-0210ISQ 7000 with VPIISQ 7000 Front Cover (with VPI)IR120191-0200ISQ 7000 no VPIISQ 7000 Front Cover (with VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0200TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0200TSQ 9000 Fron tOVEITSQ 9000 Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 Front VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 no VPIFront	Right Side Panel-ISQ 7000	1R120191-1020	ISQ 7000	
Top Manifold Cover (glass)IR120401-3000ISQ AllISQ 7610 Front Cover (with VPI)IR120191-0110ISQ 7610 with VPIISQ 7610 Front Cover (no VPI)IR120191-0210ISQ 7610 no VPIISQ Front Cover (with VPI)IR120407-0100ISQ DISQ Front Cover (with VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120189-0121TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 mov HVITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITS	Left Side Panel-ISQ	1R120411-0001	ISQ All	
ISQ 7610 Front Cover (with VPI)IR120191-0110ISQ 7610 with VPIISQ 7610 Front Cover (no VPI)IR120191-0210ISQ 7610 no VPIISQ Front Cover (with VPI)IR120407-0100ISQ DISQ Front Cover (with VPI)IR120191-0100ISQ QDISQ 7000 Front Cover (no VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0200TSQ AllLeft Hand Sub Panel-TSQIR120541-0200TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9610 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Evo Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Fr	Top Cover Panel-ISQ	1R120412-0001	ISQ All	
ISQ 7610 Front Cover (no VPI)IR120191-0210ISQ 7610 no VPIISQ Front Cover (with VPI)IR120407-0100ISQ, ISQ LTISQ Front Cover (No VPI)IR120407-0200ISQ QDISQ 7000 Front Cover (with VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Class Cover-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 no VPITSQ 9610 Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120548-0300TSQ 8000 Evo mo VPITSQ 9000 Front Cover (with VPI)IR120548-0300TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120548-0200TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 mo VPITSQ 9000 Front Cover (no VPI)IR120548-0200TSQ 9000 mo VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door Hinge SupportIR1	Top Manifold Cover (glass)	1R120401-3000	ISQ All	-
ISQ Front Cover (with VPI)1R120407-0100ISQ, ISQ LTISQ Front Cover (No VPI)1R120407-0200ISQ QDISQ 7000 Front Cover (with VPI)1R120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)1R120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQ1R120542-3000TSQ AllSmall Top Glass Cover-TSQ1R120541-0200TSQ AllLeft Hand Front Panel-TSQ1R120541-0003TSQ AllLeft Hand Sub Panel-TSQ1R120541-0003TSQ AllTSQ 9610 Front Cover (with VPI)1R120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)1R120192-0120TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)1R120192-0120TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 9000 Evo Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Sup	ISQ 7610 Front Cover (with VPI)	1R120191-0110	ISQ 7610 with VPI	
ISQ Front Cover (No VPI)IR120407-0200ISQ QDISQ 7000 Front Cover (with VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllRight Side Panel-TSQIR120541-0003TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120192-0220TSQ 9610 with NVAEITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 ith VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR12045-0001AllFro	ISQ 7610 Front Cover (no VPI)	1R120191-0210	ISQ 7610 no VPI	
ISQ 7000 Front Cover (with VPI)IR120191-0100ISQ 7000 with VPIISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120541-0200TSQ AllRight Side Panel-TSQIR120541-0003TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllTop Cover Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120192-0120TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Evo Front Cover (with VPI)IR120192-0100TSQ 9000 Evo with VPITSQ 9000 Evo Front Cover (no VPI)IR120192-0100TSQ 9000 Evo with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 Evo mot VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 Evo mith VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 mot VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 mot VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 mot VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 mot VPITSQ 9000 Front Cover (no VPI)IR120145-0001AllFront Door HingeIR3434-5000AllFront Door Hinge SupportIR120445-0001AllFront Door LatchIR76483-3000Al	ISQ Front Cover (with VPI)	1R120407-0100	ISQ, ISQ LT	
ISQ 7000 Front Cover (no VPI)IR120191-0200ISQ 7000 no VPILarge Top Glass Cover-TSQIR120542-3000TSQ AllSmall Top Glass Cover-TSQIR120401-3000TSQ AllRight Side Panel-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0001TSQ AllTop Cover Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with AEI VPI)IR120192-0120TSQ 9610 no VPITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 Evo no VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door Hinge SupportIR120445-0001AllFront Door LatchIR76483-3000AllRight Side Panel FeetIR3666-0207All	ISQ Front Cover (No VPI)	1R120407-0200	ISQ QD	-
Large Top Glass Cover-TSQ1R120542-3000TSQ AllSmall Top Glass Cover-TSQ1R120401-3000TSQ AllRight Side Panel-TSQ1R120541-0200TSQ AllLeft Hand Front Panel-TSQ1R120541-0003TSQ AllLeft Hand Sub Panel-TSQ1R120541-0001TSQ AllLeft Hand Sub Panel-TSQ1R120541-0001TSQ AllTop Cover Panel-TSQ1R120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)1R120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with AEI VPI)1R120192-0220TSQ 9610 no VPITSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (with VPI)1R120548-0400TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0100TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllRight Side Panel Feet1R3666-0207All	ISQ 7000 Front Cover (with VPI)	1R120191-0100	ISQ 7000 with VPI	-
Small Top Glass Cover-TSQIR120401-3000TSQ AllRight Side Panel-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0001TSQ AllTop Cover Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with VPI)IR120192-0220TSQ 9610 with NVAEITSQ 9610 Front Cover (with AEI VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo with VPITSQ 8000 Evo Front Cover (no VPI)IR120548-0200TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 mo VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door LatchIR76483-3000AllRight Side Panel FeetIR3666-0207All	ISQ 7000 Front Cover (no VPI)	1R120191-0200	ISQ 7000 no VPI	-
Right Side Panel-TSQIR120541-0200TSQ AllLeft Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0001TSQ AllTop Cover Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (with AEI VPI)IR120192-0200TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120192-0120TSQ 9610 with NVAEITSQ 9610 Front Cover (with AEI VPI)IR120189-0121TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (with VPI)IR120548-0400TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door Hinge SupportIR120445-0001AllFront Door LatchIR76483-3000AllRight Side Panel FeetIR3666-0207All	Large Top Glass Cover-TSQ	1R120542-3000	TSQ All	-
Left Hand Front Panel-TSQIR120541-0003TSQ AllLeft Hand Sub Panel-TSQIR120541-0500TSQ AllTop Cover Panel-TSQIR120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (no VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120548-0300TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 9000 Front Cover (with VPI)IR120548-0200TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (with VPI)IR120192-0200TSQ 9000 no VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door Hinge SupportIR120445-0001AllFront Door LatchIR76483-3000AllRight Side Panel FeetIR3666-0207All	Small Top Glass Cover-TSQ	1R120401-3000	TSQ All	-
Left Hand Sub Panel-TSQ1R120541-0500TSQ AllTop Cover Panel-TSQ1R120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)1R120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (no VPI)1R120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)1R120192-0120TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 Evo with VPITSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	Right Side Panel-TSQ	1R120541-0200	TSQ All	-
Left Hand Sub Panel-TSQ1R120541-0500TSQ AllTop Cover Panel-TSQ1R120541-0001TSQ AllTSQ 9610 Front Cover (with VPI)1R120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (no VPI)1R120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)1R120189-0121TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 Evo with VPITSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	Left Hand Front Panel-TSQ	1R120541-0003	TSQ All	
TSQ 9610 Front Cover (with VPI)IR120192-0120TSQ 9610 with VPITSQ 9610 Front Cover (no VPI)IR120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)IR120189-0121TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)IR120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)IR120548-0400TSQ 8000 Evo with VPITSQ Duo Front Cover (no VPI)IR120548-0200TSQ 8000 Evo with VPITSQ 9000 Front Cover (with VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)IR120192-0200TSQ 9000 no VPIFront Door HingeIR3434-5000AllFront Door LatchIR76483-3000AllRight Side Panel FeetIR3666-0207All	Left Hand Sub Panel-TSQ	1R120541-0500	TSQ All	- Cover Components
TSQ 9610 Front Cover (no VPI)1R120192-0220TSQ 9610 no VPITSQ 9610 Front Cover (with AEI VPI)1R120189-0121TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ Duo Front Cover (with VPI)1R120548-0200TSQ DuoTSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	Top Cover Panel-TSQ	1R120541-0001	TSQ All	-
TSQ 9610 Front Cover (with AEI VPI)1R120189-0121TSQ 9610 with NVAEITSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ Duo Front Cover1R120548-0200TSQ DuoTSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 9610 Front Cover (with VPI)	1R120192-0120	TSQ 9610 with VPI	-
TSQ 8000 Evo Front Cover (with VPI)1R120548-0300TSQ 8000 Evo no VPITSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ Duo Front Cover1R120548-0200TSQ DuoTSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 9610 Front Cover (no VPI)	1R120192-0220	TSQ 9610 no VPI	-
TSQ 8000 Evo Front Cover (no VPI)1R120548-0400TSQ 8000 Evo with VPITSQ Duo Front Cover1R120548-0200TSQ DuoTSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 9610 Front Cover (with AEI VPI)	1R120189-0121	TSQ 9610 with NVAEI	-
TSQ Duo Front Cover1R120548-0200TSQ DuoTSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 8000 Evo Front Cover (with VPI)	1R120548-0300	TSQ 8000 Evo no VPI	-
TSQ 9000 Front Cover (with VPI)1R120192-0100TSQ 9000 with VPITSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 8000 Evo Front Cover (no VPI)	1R120548-0400	TSQ 8000 Evo with VPI	-
TSQ 9000 Front Cover (no VPI)1R120192-0200TSQ 9000 no VPIFront Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ Duo Front Cover	1R120548-0200	TSQ Duo	
Front Door Hinge1R3434-5000AllFront Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 9000 Front Cover (with VPI)	1R120192-0100	TSQ 9000 with VPI	
Front Door Hinge Support1R120445-0001AllFront Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	TSQ 9000 Front Cover (no VPI)	1R120192-0200	TSQ 9000 no VPI	
Front Door Latch1R76483-3000AllRight Side Panel Feet1R3666-0207All	Front Door Hinge	1R3434-5000	All	
Right Side Panel Feet1R3666-0207All	Front Door Hinge Support	1R120445-0001	All	
	Front Door Latch	1R76483-3000	All	
Chassis Foot 1R3666-0206 All	Right Side Panel Feet	1R3666-0207	All	1
	Chassis Foot	1R3666-0206	All	1

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Small Turbomolecular Pump	1R119268-0004	ISQ All	
Medium Turbomolecular Pump	1R119268-0003	TSQ 9000	
Large Turbomolecular Pump	1R119268-0002	All	
Rough Pump (RV3)	1R76505-3007	All	Dump Components
Dry Scroll Pump	120256-0001	All	— Pump Components
Vacuum Pump Oil	A0301-15101	All	
Convectron Gauge	A0105-00501	All	1
Foreline Adapter	1R119244-0025	All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
Bushing in the Source Exchange	1R120406-2203	All	
Seal in the Source Exchange Tool	1R120406-2204	All	
Clip in the Source Exchange Tool	1R76483-2102	All	
Source Exchange Tool	1R120406-2000	All	
Source Holder	1R120471-0001	All except AEI MS Systems	
Source Holder, NVAEI MS systems	1R120189-2023	TSQ 9610 with NVAEI	
T30 Torxhead Key	1R3812-5T30	All	
MS Toolkit, which includes:	1R120467-0001	All No VPI MS Systems	
Source Removal Tool (Small)	1R120406-2250	All No VPI MS Systems	
Column Measurement Tool	1R120461-0010	All	
T10 Torxhead Key	1R3812-5T10	All	
T20 Torxhead Key	1R3812-5T20	All	
• Forceps, 8 in.	1R76360-0008	All	MS System Tools
• Wrench, Open-Ended, 1/4-in., 5/16-in.	1R76360-0109	All	
• Wrench, Open-Ended, 3/8 in., 7/16-in.	1R76360-0108	All	
VPI MS Toolkit, which includes:	1R120467-0002	ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo	
Source Removal Tool (Small)	1R120406-2250	ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo	
Column Measurement Tool	1R120461-0010	All	
Source Exchange Tool	1R120406-2000	ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo	
Source Holder	1R120471-0001	All except NVAEI MS Systems	1
• T10 Torxhead Key	1R3812-5T10	All	1
• T20 Torxhead Key	1R3812-5T20	All	1
• T30 Torxhead Key	1R3812-5T30	All]

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
• Forceps, 8 in.	1R76360-0008	All	
• Wrench, Open-Ended, 1/4-in., 5/16-in.	1R76360-0109	All	
• Wrench, Open-Ended, 3/8 in., 7/16-in.	1R76360-0108	All	
ISQ 7000, ISQ 7610, TSQ 9000, and TSQ 9610 VPI MS Toolkit, includes everything in VPI MS Toolkit and:	1R120467-0004	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610 except TSQ 9610 with NVAEI	
Source Plug	1R120589-2000	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610 except TSQ 9610 with NVAEI	
Source Plug O-ring	1R3816-202	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610 except TSQ 9610 with NVAEI	
Source Plug Holder	1R120589-1050	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610 except TSQ 9610 with NVAEI	MS System Tools
NVAEI MS Toolkit	1R120189-6200	TSQ 9610 with NVAEI	
NVAEI Source Exchange Tool	1R120189-0003	TSQ 9610 with NVAEI	
NVAEI Source Holder	1R120189-2303	TSQ 9610 with NVAEI	
NVAEI Manual Source Removal Tool	1R120189-1017	TSQ 9610 with NVAEI	
NVAEI Source Plug	1R120189-6011	TSQ 9610 with NVAEI	
Column Measurement Tool	1R120461-0010	All	
T10 Torxhead Key	1R3812-5T10	All	
• T20 Torxhead Key	1R3812-5T20	All	
• Forceps, 8 in.	1R76360-0008	All	
• Wrench, Open-Ended, 1/4-in., 5/16-in.	1R76360-0109	All	
• Wrench, Open-Ended, 3/8 in., 7/16-in.	1R76360-0108	All	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
CI Reagent Gas Flow Module	1R23331-0092	All	
CI Ion Volume	1R120404-4112	All	
CI Ion Source Cartridge Assembly (Low Activity)	1R120404-4500	All	
CI Gas Line	1R120438-0020	All	
EI/CI Ion Volume	1R120404-4113	All	
Hydrogen Ion Volume	1R120404-4115	All	
Dynode and Multiplier Power Supply and Cables (EI/CI)	1R120361-0007	All	_
Ion Gauge Option	1R120560-0020	All	
Ion Gauge Controller PCB	1R120485-0120	All	
Ion Gauge Mount	1R120416-0002	All	
Ion Gauge Tube Shield	1R119605-0012	All	
Extended Capacity Turbomolecular Pump	1R119268-0002	All	
Dust Filters	1R120442-1000	All	
Direct Insert Probe (without case)	1R120406-4000	All VPI Systems	- Upgrade Equipment
Direct Exposure Probe (without case)	1R120406-5000	All VPI Systems	
Direct Insert Probe Controller Kit	1R119300-0500	All VPI System	
Direct Exposure Probe Controller Kit	1R119300-0600	All VPI System	
AEI Ion Source Cartridge, which includes:	1R120602-1300	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
AEI Ion Volume	1R120602-0104	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
• Dual Filament	1R120404-1940	All	
AEI Filament Spacer	1R120602-0103	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
• AEI Filament Retainer	1R120602-0102	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
• AEI Ion Volume-Repeller Insulator	1R120602-0106	IISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
• AEI Filament Shield	1R120602-0105	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
AEI Ion Source Sleeve	1R120602-0101	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
• AEI Filament Insulator	1R120602-0108	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
Lens 3/RF Lens	1R120574-0103	All	
• Lens 2	1R120404-1140	All	
NeverVent AEI Ion Source Cartridge, which includes:	1R120189-6030	TSQ 9610 with NVAEI	
NeverVent AEI Ion Volume	1R120189-2001	TSQ 9610 with NVAEI	
• NeverVent AEI Dual Filament Assembly	1R120189-6000	TSQ 9610 with NVAEI	
NeverVent AEI Ion Cartridge Sleeve	1R120189-2000	TSQ 9610 with NVAEI	
• NeverVent AEI Filament Insulator	1R120189-2036	TSQ 9610 with NVAEI	
NeverVent AEI Cap	1R120189-2035	TSQ 9610 with NVAEI	
NeverVent AEI Filament Spring	1R120189-3002	TSQ 9610 with NVAEI	
• AEI Filament Shield	1R120602-0105	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
Lens 3/RF Lens	1R120574-0103	All	
• Lens 2	1R120404-1140	All	
AEI Front Plate Assembly	1R120191-0014	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
AEI Case Insert	1R120602-1100	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	Upgrade Equipment
AEI Case	1R76524-0020	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
AEI Magnet Retainer Bracket	1R120602-1104	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
AEI Magnet	1R76253-2000	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	
AEI Source Magnet Assembly	1R120602-0200	ISQ 7000, ISQ 7610, TSQ 9000, TSQ 9610	

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
MS Installation Kit	1R120459-0001	All	
Test Column	26070-1300	All	
ISQ Specification Standards	1R120150-PERF	ISQ All	
ISQ Series Instrument Control Software	1R120479-0003	ISQ All	
TSQ Specification Standards	1R120540-PERF	TSQ All	
TSQ Series Software	1R120579-0031	TSQ All	
Fuse Kit, which includes:	1R120458-0001	All	
 Fuse, Fast Acting 500 mA, 250 V, 5 × 20 (pkg of 2) 	-		
• Fuse, Fast Acting 1.0 A, 250 V, 5 × 20 (pkg of 2)	-		Other
• Fuse, Fast Acting 1.6 A, 250 V, 5 × 20 (pkg of 2)	-		
• Fuse, Fast Acting 2.5 A, 250 V, 5 × 20 (pkg of 2)	-		
• Fuse, Fast Acting 5.0 A, 250 V, 5 × 20 (pkg of 2)	-		
• Fuse, Fast Acting 6.3 A, 250 V, 5 × 20 (pkg of 2)	-		
Consumables Kit	1R120468-CONSU M	All	
Aluminum Oxide	1R32000-60340	All	
Cardinal Health CP100 Nitrile Cleanroom Gloves (Recommended to avoid contaminating instrument)	Fisher Scientific Catalog Number	All	
Size X-Small (Case)	19-120-2947	All	
Size Small (Case)	19-120-2947A	All	Gloves
Size Medium (Case)	19-120-2947B		
Size Large (Case)	19-120-2947C		
Size X-Large (Case)	19-120-2947D		

Component	Thermo Scientific Part Number	Used on Instrument(s)	Category
MS Screw Kit, which includes:	1R120480-0014	All	
• M3 x 6 mm Screw (pkg of 10)			
• M3 x 8 mm Screw (pkg of 10)			
• M3 x 10 mm Screw (pkg of 5)			
• M3 x 12 mm Screw (pkg of 5)			
• M3 x 16 mm Screw (pkg of 5)			
• M4 x 16 mm Screw (pkg of 5)			
• M4 x 20 mm Screw (pkg of 10)			
• M4 x 10 mm Screw (pkg of 10)			
• M3 x 8 mm Screw for Hinges (pkg of 5)			
• M4 x 8 mm Screw for Power Supply (pkg of 5)	-		
• M3 Locking Nut for Front Door Latch (pkg of 2)			
• M3 x 25mm Standoff for RF Board (pkg of 2)			

* All of the screws in this table must be cleaned before you install them inside the vacuum manifold. See *Cleaning Durable Components* in the *ISQ 7000 Hardware Manual* or *TSQ 9000 Hardware Manual* for details.

Calibration Gas Components

You can purchase the following calibration gas components for the ISQ or TSQ GS-MS instruments mass spectrometer. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for information about installing these components or Table 1 for a list of available components.

The following calibration gas components can be replaced on the ISQ or TSQ GS-MS instrument.

- Transfer Line Assembly (New)
- Transfer Line (Previous Model)
- Single-Level Calibration Gas Module
- Dual-Level Calibration Gas Module
- Calibrant Reservoir Kit
- Gas Mixing Chamber (Legacy)
- Gas Mixing Chamber
- Source Gas Tube
- Calibration Gas Line
- FC-43 Calibration Compound
- 1 fg/µL OFN IDL Kit
- 5 fg/µL OFN in Iso-octane






Figure 1. Replaceable Calibration Gas Components-ISQ, ISQ QD, IST LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo





Transfer Line Assembly (New)


Quantity: Each

System: ISQ 7000,TSQ 9000, ISQ 7610, TSQ 9610





Thermo Scientific Part Number

1R120191-0020



NOTE: To replace all the M4 x 10 mm screws on the transfer line, order two of PN 1R76913-0410. To replace the o-ring on the transfer line, order PN 1R3814-223. PN 1R120402-0001 includes one o-ring initially.





NOTE: To replace all the M4 x 10 mm screws on the calibration controller, order two of PN 1R76913-0410.

Calibrant Reservoir Kit	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120433-0001
Gas Mixing Chamber (Legacy)	Quantity: Each
	System: ISQ, TSQ 8000

Thermo Scientific Part Number

1R120404-1201

NOTE: To replace all the M3 x 8 mm screws on the gas mixing chamber, order two of PN 1R76913-0306. This part is replaced with 1R120574-0054 and 1R120574-0055.

Gas Mixing Chamber	Quantity: Each
	System: ISQ, ISQ LT, ISQ QD, TSQ 8000, TSQ 8000 Evo, TSQ Duo

Thermo Scientific Part Number

1R120574-0054

NOTE: To replace all the M3 x 8 mm screws on the gas mixing chamber, order two of PN 1R76913-0308.

Note The Gas Mixing Chamber (PN 1R120574-0054), two Source Gas Tubes (PN 1R120574-0055), and two M3 × 8 mm screws (PN 1R176913-0308) can be ordered as a kit. To purchase the Gas Mixing Chamber Kit, order PN 1R120574-KIT.

Source Gas Tube	Quantity: Each
	System: ISQ, ISQ LT,
	ISQ QD, TSQ 8000,
0	TSQ 8000 Evo, TSQ Duo
Thermo Scientific Part Number	1R120574-0055
Calibration Gas Line	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120438-0010
FC-43 Calibration Compound	Quantity: Each
TREAD ID-Story	System: All
Thermo Scientific Part Number	1R50010-30059
1 fg/µL OFN IDL Kit	Quantity: Each
	System: ISQ 7000,TSQ 9000, ISQ 7610, TSQ 9610
Thermo Scientific Part Number	1R76310-0101K


5 fg/µL OFN in Iso-octane	Quantity: Each
	System: ISQ 7000,TSQ 9000, ISQ 7610, TSQ 9610
Thermo Scientific Part Number	1R76310-0105

Collision Gas Components

You can purchase the following collision gas components for the TSQ GC-MS systems. Be sure to reference the component's part number when placing an order with your local Sales/Service Representative. See the *TSQ 9000 Hardware Manual* for information about installing these components or Table 1 for a comprehensive list of available components.


The following collision gas components can be replaced on your mass spectrometer.

- Collision Gas Module
- Tube Fitting, Collision Gas Module to Manifold
- Collision Gas Tube
- Tee, 1/8 in.



NOTE: This PN includes the module cable (PN 1R119378-0071), tubing (PN 1R120557-0030), and 1/8 in. tee (PN 1R76256-0126).

Tube Fitting, Collision Gas Module to Manifold	Quantity: Each
	TSQ All
Thermo Scientific Part Number	1R76256-0045



Column Components

You can purchase the following column components for the ISQ or TSQ GC-MS instrument. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manuals for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following column components can be replaced on the ISQ or TSQ GC-MS instrument.

- Graphite Vespel Ferrule for 0.18 mm Columns
- Graphite Vespel Ferrule for 0.25 mm Columns
- Graphite Vespel Ferrule for 0.32 mm Columns
- Graphite Vespel Ferrule for 0.53 mm Columns
- 2-hole Graphite Vespel Ferrule for <0.32mm Columns
- No-Hole Graphite Vespel Ferrule
- SilTite Ferrule for 0.25 mm Columns
- SilTite Ferrule for 0.32mm Columns
- SilTite Ferrule for 0.53mm Columns
- Nut for SilTite Ferrules
- Nickel-Coated Nut for Graphite Vespel Ferrules
- Spring-loaded Transfer Line Nut
- Graphite Ferrules for 0.1–0.25 mm Columns used with Spring-loaded Transfer Line Nut

Note The parts in this chapter are compatible with all GC-MS systems.

Graphite Vespel Ferrule for 0.18 mm Columns




Thermo Scientific Part Number

1R76458-2016

Quantity: Pkg of 10

Graphite Vespel Ferrule for 0.25 mm Columns	Quantity: Pkg of 10
Thermo Scientific Part Number	29033496
Graphite Vespel Ferrule for 0.32 mm Columns	Quantity: Pkg of 10
Thermo Scientific Part Number	1R76458-2019
Graphite Vespel Ferrule for 0.53 mm Columns	Quantity: Pkg of 10
Thermo Scientific Part Number	1R76458-2020
2-hole Graphite Vespel Ferrule for <0.32mm Columns	Quantity: Pkg of 10
AII	
Thermo Scientific Part Number	1R76458-2018



SilTite Ferrule for 0.53mm Columns	Quantity: Pkg of 10
Thermo Scientific Part Number	1R76458-2026
Nut for SilTite Ferrules	Quantity: Pkg of 5
Thermo Scientific Part Number	1R76458-2001
Nickel-Coated Nut for Graphite Vespel Ferrules	Quantity: Pkg of 5
Thermo Scientific Part Number	1R76256-0060

Spring-loaded Transfer Line Nut	Quantity: Each
Thermo Scientific Part Number	1R120434-0010
Graphite Ferrules for 0.1–0.25 mm Columns used with Spring-loaded Transfer Line Nut	Quantity: Pkg of 10
Thermo Scientific Part Number	290VA191

Fan and Filter Components

You can purchase the following filter components for the ISQ and TSQ GC-MS systems. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following fan and filter components can be replaced.

- Intake Filters/RF Shields
- Chassis Cooling Fan
- Oil Mist Filter
- Turbo Cooling Fan
- Fan Plenum
- Left Rear Wire Mesh Filter
- Front Wire Mesh Filter
- Air Deflector

Intake Filters/RF Shields

Quantity: Each


System: All



Thermo Scientific Part Number

1R76475-5000

NOTE: The ISQ Series instrument has two intake filters, so to replace them, order two of PN 1R76475-5000. To replace all the M3 x 8 mm screws on the filters, order eight of PN 1R76913-0308.



Thermo Scientific Part Number

NOTE: To replace all of the nylon rivets on the chassis cooling fan, order four of PN 1R3326-5000.


Oil Mist Filter


Quantity: Each


System: All

Thermo Scientific Part Number

1R76505-0036








Ion Source Components

The following ion source components can be replaced on the ISQ and TSQ GC-MS systems. Be sure to include the component's part number when placing your order. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

- Ion Source Block (Old)
- ISQ 7000, ISQ 7610, TSQ 9000 and TSQ 9610 Ion Source Block
- TSQ 9610 with NVAEI Ion Source Block
- Repeller Plate
- Source-Repeller Thumbscrew (Old)
- ISQ 7000, ISQ 7610, TSQ 90000 and TSQ 9610 Source-Repeller Thumbscrew
- Thumbscrew-Repeller Spacer
- Insulating Spacer
- Lens/Source Heater
- Lens Heater Block
- Lens Heater Block Grounding Strap
- Lens Plate and Springs
- Source to Ion Guide Spacer
- EI Ion Source Cartridge (Low Activity)
 - Ion Cartridge Sleeve
 - EI Ion Volume (Low Activity)
 - Ion Volume-Repeller Insulator
 - Repeller (Low Activity)
 - Ion Volume Locking Ring
 - Repeller Spring
 - Repeller Nut
 - Lens 3/RF Lens
 - Lens 1
 - Lens 2





Ion Source Block (Old)	Quantity: Each
	System: ISQ, ISQ LT, ISQ QD, TSQ Duo, TSQ 8000, TSQ 8000 Evo
Thermo Scientific Part Number	1R120404-1311
ISQ 7000, ISQ 7610, TSQ 9000 and TSQ 9610 Ion Source Block	Quantity: Each
	System: ISQ 7000, ISQ 7610, TSQ 9000, and TSQ 9610 without AEI
Thermo Scientific Part Number	1R120191-2000
TSQ 9610 with NVAEI Ion Source Block	Quantity: Each
TSO 9610 with NVAEI Ion Source Block	Quantity: Each System: TSQ 9610 with NVAEI

Repeller Plate	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1312
Source-Repeller Thumbscrew (Old)	Quantity: Each
Source-neperier rinumbscrew (Olu)	Quantity: Each
	System: ISQ
Thermo Scientific Part Number	1R120404-1400

Thermo Scientific Part Number	1R120404-1400	
Source-Repeller Thumbscrew (Old)	Quantity: Each	
	System: ISQ, ISQ LT, ISQ QD, TSQ 8000, TSQ 8000 Evo, TSQ Duo	

Thermo Scientific Part Number	1R120564-0004
Note If your instrument has a source insertion guide, orde	er two of PN 1R120564-0003. If your instrument was not
manufactured with a source insertion guide, order two of	PN 1R120404-1400. See Manifold Components to
determine if your instrument was manufactured with a so	urce insertion guide.

ISQ 7000, ISQ 7610, TSQ 90000 and TSQ 9610 Source-Repeller Thumbscrew	Quantity: Each
Æ	System: ISQ 7000, ISQ 7610, TSQ 9000, and TSQ 9610
Thermo Scientific Part Number	1R120602-0206

Thumbscrew-Repeller Spacer	Quantity: Each
	System: All
0	

Thermo Sci	entific Part Numb	er		1RTSLI053A	

NOTE: There are two thumbscrew-repeller spacers on the ion source, so to replace them, order two of PN 1RTSLI053A.


Insulating Spacer	Quantity: Each
	System: All except TSQ 9610 with AEI


Thermo Scie	ntific Part Number	1RTSLI050A

NOTE: There are two insulating spacers on the ion source, so to replace them, order two of PN 1RTSLI050A.

Lens/Source HeaterQuantity: EachSystem: AllSystem: AllThermo Scientific Part NumberIR120404-1330

NOTE: There are two lens/source heaters on the ion source, so to replace them, order two of PN 1R120404-1330. To replace the M4 x 20 mm screws on each lens/source heater, order two of PN 1R76913-0420.





NOTE: To replace all the M3 x 6 mm screws on the grounding strap, order two of PN 1R76913-0305.

Lens Plate and Springs	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1750

NOTE: To replace the M3 x 6 mm screw on the lens plate and springs, order PN 1R76913-0306. To replace the copper washer, order PN 1R76483-5003.

Source to Ion Guide Spacer	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1340

Ion Source Cartridge Components

The ion source cartridge consists of the following components. See the hardware manual for information about installing these components or Table 1 for a comprehensive list of available components.

- EI Ion Source Cartridge (Low Activity)
- Ion Cartridge Sleeve
- EI Ion Volume (Low Activity)
- Ion Volume-Repeller Insulator
- Repeller (Low Activity)
- Ion Volume Locking Ring
- Repeller Spring
- Repeller Nut
- Lens 3/RF Lens
- Lens 1
- Lens 2

El Ion Source Cartridge (Low Activity)	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-4100
Ion Cartridge Sleeve	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1105

El Ion Volume (Low Activity)	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-4111
Ion Volume-Repeller Insulator	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1114
Repeller (Low Activity)	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1161
las Valuma Lasking Ding	Quantity Fach
Ion Volume Locking Ring	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1118

Repeller Spring	Quantity: Pkg of 5
	System: All
Thermo Scientific Part Number	1R76485-1000K
NOTE: Although the ion source only contains one repeller spring, you v	will receive a package of five when you order PN 1R76485-1000K.
Repeller Nut	Quantity: Each
0	System: All
Thermo Scientific Part Number	1R120404-1120
Lens 2	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-1140
Lens 3/RF Lens	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120574-0103

NOTE: This part replaces PN 1R120404-1150.

	System: All
Thermo Scientific Part Number	1R120574-0139

NOTE: This part replaces PN 1R120404-1130.


Electron Multiplier Components

You can purchase the following electron multiplier components for the ISQ and TSQ GC-MS system. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following electron multiplier components can be replaced.

- Electron Multiplier (Old)
- Electron Multiplier Plate
- Electron Multiplier, ISQ 7610 and TSQ 9610
- Dynode Feedthrough
- Anode Feedthrough Assembly
- Detector Assembly (Dynode, Electron Multiplier, and Lead Set)


Table 3. Replaceable Components of the Electron Multiplier



Electron Multiplier (Old)

Quantity: Each

System: ISQ, ISQ QD, ISQ LT, ISQ 7000, TSQ 8000, TSQ 8000 Evo, TSQ Duo, TSQ 9000




Thermo Scientific Part Number

1R120405-1000

NOTE: To replace the M4 x 10 mm screws, which secure the electron multiplier plate to the manifold, order four of PN 1R76913-0306.

Electron Multiplier, ISQ 7610 and TSQ 9610	Quantity: Each
	System: ISQ7610 and TSQ 9610
Thermo Scientific Part Number	1R76022-14634
Dura da Fasaldaran k	
Dynode Feedthrough	Quantity: Each
Large Feedthrough	System: All
0-Ring	
Thermo Scientific Part Number	1R120610-0040

NOTE: There are two large feedthroughs on the ISQ Series instrument, so to replace them, order two of PN 1R120610-0040. To replace all the M4 x 10 mm screws on each large feedthrough, order three of PN 1R76913-0410. To replace the o-ring, order PN 1R3814-113.



NOTE: To replace all the M3 x 6 mm screws on the anode feedthrough, order two of PN 1R76913-0306. To replace the o-ring, order PN 1R3814-110. To replace the standoff, order PN 1R77005-3010. To replace the anode feedthrough, order PN 1R120610-0050.


Detector Assembly (Dynode, Electron Multiplier, and Lead Set)	Quantity: Each	
	System: All	
Thermo Scientific Part Number	1R120379-0001	


Dual Filament Components

You can purchase the following dual filament components for the ISQ and TSQ GC-MS systems. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following dual filament components can be replaced.

- Dual Filament
- Filament Retaining Spring and Spacer
- Figure 3. Replaceable Components of the Dual Filament





NOTE: PN 1R120404-1405 contains the filament retaining spring, spacer, and screw. To replace only the M3 x 6 mm screw on the filament retaining spring, order PN 1R76913-0306.

Analyzer Components

You can purchase the following analyzer components for the ISQ or TSQ GC-MS instruments. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following analyzer components can be replaced on the ISQ or TSQ GC-MS instrument.

- ISQ EI Analyzer Tray with Quad Assembled (ISQ Classic, ISQ QD, ISQ LT)
- ISQ Analyzer Tray (Old)
- ISQ 7000 and ISQ 7610 Analyzer Tray Assembled
- ISQ Analyzer Tray (New)
- Ion Guide
- Ion Guide Clamp
- Q1 Assembly with Wires and Tested
- ISQ Exit Endcap
- Entrance Endcap
- Q1 Quad Entrance Lens
- Ion Guide Wire Set
- Collision Cell Assembly
- Q3 Quad Assembly
- Q3 Exit Lens Assembly
- Quad 1 Wire Set
- Quad 3 Wire Set
- TSQ Analyzer Tray with Source Block and Ion Guide (Old)
- TSQ 9610 with NVAEI Analyzer Tray Assembled
- TSQ 9610 with NVAEI Analyzer Tray Assembled
- Wire Routing Clip
- Q1 Exit Lens Wire
- Q3 Entrance Lens Wire
- Q3 Exit Lens Wire






Figure 4. Replaceable Components of the ISQ GC-MS Instrument Analyzer

* See Board Components for information about ordering the source interface board. ** See for information about ordering the ion source.





* See Ion Source Components for information about ordering the ion source.

SQ EI Analyzer Tray with Quad Assembled (ISQ Classic, ISQ QD, ISQ LT)	Quantity: Each
	System: ISQ, ISQ LT, ISQ QD
hermo Scientific Part Number	1R120404-0003
SQ Analyzer Tray (Old)	Quantity: Each
	System: ISQ, ISQ LT, ISQ QD
Thermo Scientific Part Number	1R120404-2000
	1R120404-2000 Quantity: Each
Thermo Scientific Part Number ISQ 7000 and ISQ 7610 Analyzer Tray Assembled	Quantity: Each System: ISQ 7000 and

SQ Analyzer Tray (New)	Quantity: Each
	System: All ISQ systems
Thermo Scientific Part Number	1R120404-2000
on Guide	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120404-3100



NOTE: To replace all the M3 x 16 mm screws on the ion guide clamp, order two of PN 1R76913-0316.


Q1 Assembly with Wires and Tested	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120542-2560

ISQ Exit Endcap	Quantity: Each
	System: ISQ All
Thermo Scientific Part Number	1R120404-3230
OTE: To replace all the M4 x 6 mm screws on the exit endcap,	order two of PN 1R76913-0406.
Entrance Endcap	Quantity: Each
	System: All


Thermo Scientific Part Number

1R120404-3211

NOTE: To replace all the M4 x 6 mm screws on the entrance endcap, order two of PN 1R76913-0406.




NOTE: To replace all the M3 x 16 mm screws on the quad entrance lens, order two of PN 1R76913-0316.




NOTE: To replace the M3 nut on the ion guide wire set, order PN 1R76904-0003.

Collision Cell Assembly Quantity: Each System: TSQ All **Thermo Scientific Part Number** 1R120574-0161 NOTE: To replace all the M4 x 6 mm screws to mount the collision cell to the analyzer tray, order four of PN 1R76913-0306. **Q3 Quad Assembly** Quantity: Each System: TSQ All **Thermo Scientific Part Number** 1R120542-2570



NOTE: To replace all the M4 x 6 mm screws on the Q3 exit lens assembly to mount it to the tray, order two of PN 1R76913-0406.



NOTE: To replace all the M3 x 5 mm screws on the quad wire set, order four of PN 1R76913-0305. To replace all the copper washers, order eight of PN 1R76483-5003.

Quad 3 Wire Set

Quantity: Pkg of 2

System: TSQ All



Thermo Scientific Part Number

1R120480-0006


NOTE: To replace all the M3 x 5 mm screws on the quad wire set, order four of PN 1R76913-0305. To replace all the copper washers, order eight of PN 1R76483-5003.

TSQ Analyzer Tray with Source Block and Ion Guide (Old)	Quantity: Each
	Systems: TSQ 8000, TSQ 8000 Evo, TSQ Duo

1R120542-2535

SQ 9000 Analyzer Tray Assembled	Quantity: Each
	System: TSQ 9000
ermo Scientific Part Number	1R120192-0011T
ire Routing Clip	Quantity: Each
	System: TSQ All

TSQ 9610 with NVAEI Analyzer Tray Assembled	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-0011T
01 Exit Lens Wire	Quantity: Each
	System: TSQ All
	oysteni. ToQ fui



Board Components

You can purchase the following board components for the ISQ or TSQ GC-MS system. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following board components can be replaced.

- Q1 RF Board
- Q3 RF Board
- Q1/Rod Driver Board
- ISQ Lens Driver Board
- Q3 Rod Driver Board (75V)
- Q3 Rod Driver Board (75V)
- TSQ Lens Driver Board (Old)
- TSQ 9000 and TSQ 9610 Lens Driver Board
- Source Interface Board, ISQ Systems
- Source Interface Board, TSQ Systems
- Heat Shield for the AEI Source Interface Board
- Heat Shield for the Source Interface Board
- Electrometer Board
- Electrometer Shield
- Controller Interface Board
- Controller Interface Board Support Bracket
- PC Communication Board
- ISQ Distribution Board
- TSQ Distribution Board
- 20-Pin Feedthrough
- 4-Pin Feedthrough




Figure 6. Replaceable Board Components of the ISQ GC-MS System



WARNING ELECTRICAL SHOCK HAZARD: For safety reasons, the fuses on the lens driver board can only be replaced by Field Service Engineers.





Figure 7. Replaceable Board Components of TSQ GC-MS Systems



WARNING ELECTRICAL SHOCK HAZARD: For safety reasons, the fuses on the lens driver board can only be replaced by Field Service Engineers.



Q1/Rod Driver Board	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120485-0415
	1R120485-0415
NOTE: If you only need to replace the fan, order PN 1R120443-0011.	1R120485-0415 Quantity: Each
Thermo Scientific Part Number NOTE: If you only need to replace the fan, order PN 1R120443-0011. ISO Lens Driver Board	

NOTE: If you only need to replace the fan, order PN 1R120443-0002


Q3 Rod Driver Board (75V)

Quantity: Each


System: TSQ All




Thermo Scientific Part Number	1R120485-0420
TSQ Lens Driver Board (Old)	Quantity: Each
	System: TSQ 8000, TSQ 8000 Evo, TSQ Duo
Thermo Scientific Part Number	1R120485-A110




Source Interface Board, TSQ Systems	Quantity: Each
	System: TSQ All
hermo Scientific Part Number	1R120485-0210
Heat Shield for the AEI Source Interface Board	Quantity: Each
	System: TSQ 9610 with NVAEI





NOTE: To replace the M4 x 10 mm screw on the heat shield, order PN 1R76913-0410.




NOTE: To replace all the M3 x 6 mm screws on the electrometer board, order two of PN 1R76913-0306. To replace all the M3 x 16 mm screws, order two of PN 1R76913-0316.










NOTE: To replace all the M3 x 10 mm screws on the distribution board, order 26 of PN 1R76913-0310.



Thermo Scientific Part Number

NOTE: To replace all the M4 x 10 mm screws on the 20-pin feedthrough, order three of PN 1R76913-0410. To replace the o-ring, order PN 1R3814-123.

4-Pin Feedthrough Quantity: Each System: All Feedthrough Genetific Part Number

NOTE: To replace all the M4 x 10 mm screws on the 4-pin feedthrough, order three of PN 1R76913-0410. To replace the o-ring, order PN 1R3814-127.


Power Supply Components

You can purchase the following power supply components for the ISQ and TSQ GC-MS systems. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.

The following power supply components can be replaced.

- ISQ Power Supply System
- TSQ Power Supply System
- Dynode and Multiplier Power Supply and Cables (EI only)

Figure 8. Replaceable Power Supply Components of the ISQ GC-MS






Figure 9. Replaceable Power Supply Components of the TSQ GC-MS

Quantity: Each

Systems: ISQ All

NOTE: To replace all the M4 x 8 mm screws on the power supply, order two of PN 1R76937-0408.

Thermo Scientific Part Number


ISQ Power Supply System

1R120380-0001

TSQ Power Supply System

Quantity: Each


TSQ All



Thermo Scientific Part Number

1R120544-0002

NOTE: To replace all the M4 x 10 mm flat-head screws on the power supply, order PN 1R76937-0410. To replace all the M4 x 10 mm pan-head screws, order PN 1R76913-0410.



NOTE: To replace all the M4 x 10 mm screws on the dynode and multiplier power supply, order two of PN 1R76913-0410. To replace the M4 x 16 mm screw, order PN 1R76913-0416.


Vacuum Interlock Components

You can purchase the following vacuum interlock components for the any GC-MS system with a vacuum probe interlock. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manuals for your system for information about installing these components or Table 1 for a comprehensive list of available components.

The following vacuum interlock components can be replaced.

- Inner Ball Seal
- Valve Axle
- Axle Seal
- Axle Bearing
- Sealing Ball
- Outer Ball Seal
- Valve Spring
- Vacuum Interlock Knob
- Clip
- Interior O-Ring
- Exterior O-Ring
- Vacuum Interlock Microswitch Assembly
- Vacuum Interlock Assembly (Old)
- Vacuum Interlock Assembly (ISQ 7000. ISQ 7610, TSQ 9000, and TSQ 9610)
- Front Manifold Plate Assembly, TSQ 9610 with VPI AEI

Table 4. Replaceable Components of the Vacuum Interlock



Inner Ball Seal	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-1008

Valve Axle	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-1003
Axle Seal	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-2204
Axle Bearing	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-1005

Sealing Ball	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-1004
Outer Ball Seal	Quantity: Each
	System: All VPI systems
Thermo Scientific Part Number	1R120406-1002
Valve Spring	Quantity: Each
Ø	System: All VPI systems

Vacuum Interlock Knob


Quantity: Each

System: All VPI systems



Thermo Scientific Part Number	1R120406-3000


Clip		Quantity: Pkg of 5
	\mathcal{O}	System: All VPI systems
Thermo Scientific Part Number		1R76483-2102



Vacuum Interlock Assembly (Old)

Quantity: Assembly

Systems: ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo



1R120403-0001

Vacuum Interlock Assembly (ISQ 7000. ISQ 7610, TSQ 9000, and TSQ 9610)	Quantity: Assembly
	Systems: ISQ 7000. ISQ 7610, TSQ 9000, and TSQ 9610 with VPI (except TSQ 9610 with VPI AEI)


Thermo Scientific Part Number

1R120191-0012

Front Manifold Plate Assembly, TSQ 9610 with VPI AEI

Quantity: Assembly

Systems: TSQ 9610 with NVAEI



Thermo Scientific Part Number

1R120189-0012

Manifold Components

You can purchase the following manifold components for the ISQ or TSQ GC-MS system. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your system for information about installing these components or Table 1 for a comprehensive list of available components.

- Manifold Door Hinge
- Front Manifold Plate (VPI)
- Front Manifold Plate Assembly (No VPI)
- Manifold O-Ring
- Vent Valve Knob
- Source Insertion Guide
- Alignment Pin
- Front Door Alignment Pin
- Source Magnet Retainer
- Magnet Yoke
- Source Magnets
- Back Manifold Plate




Figure 10. Replaceable Components of the No VPI Manifold





Table 5. Replaceable Components of the VPI Instrument Manifold




NOTE: There are four large manifold o-rings in the ISQ Series instrument: two on the top of the vacuum manifold, one at the front and one at the back. To replace all of the manifold o-rings, order four of PN 1R3815-360.

Vent Valve Knob	Quantity: Each
	System: All
hermo Scientific Part Number	1R120403-0104
IOTE: To replace the o-ring on the vent valve knob, order PN 1R3814-110.	
Manifold Door Hinge	Quantity: Each
	System: All

Thermo Scientific Part Number	1R120565-0001
NOTE: To replace all the M4 v 10 mm percent on the manifold door bings, order three of E	NI 1076012 0410

NOTE: To replace all the M4 x 10 mm screws on the manifold door hinge, order three of PN 1R76913-0410.





NOTE: To replace all the M4 x 16 mm screws on the front manifold plate, order four of PN 1R76913-0416.

Figure 11. Replaceable Components of the No VPI Manifold

١	
Source Insertion Guide	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120564-0001
Alignment Pin	Quantity: Each
	Systems: ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo
Thermo Scientific Part Number	1R120403-0103
Front Door Alignment Pin	Quantity: Each
	Systems: ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo
Thermo Scientific Part Number	1R120564-0003

Note If your instrument has a source insertion guide, two front door alignment pins (PN 1R120564-0003) are used on the front door plate to hold the magnet yoke in place and align the door to the ion source. If your instrument was not manufactured with a source insertion guide, two alignment pins (PN 1R120403-0103) are used on the front door plate. All ISQ QD and ISQ LT instruments have 1R120564-0003 mounted on the front door plate. See Replaceable Components of the VPI Instrument Manifold to locate the source insertion guide on your instrument.


Source Magnet Retainer	Quantity: Each
	Systems: ISQ 7000; TSQ 9000

Thermo Scientific Part Number

1R120602-0205

Magnet Yoke	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120564-0002
Source Magnets	Quantity: Each
	System: All
Thermo Scientific Part Number	1R70001-98195

NOTE: There are two source magnets in the ISQ Series instrument, so to replace them, order two of PN 1R70001-98195.



Thermo Scientific Part Number

1R120403-1001

NOTE: This PN includes the plate, alignment pin (1R120403-0103) and O-ring (1R3815-360). To replace all the M4 x 16 mm screws on the back manifold plate, order four of PN 1R76913-0416.


Cover Components

You can purchase the following cover components for the ISQ and TSQ GC-MS system. Be sure to include the component's part number when placing an order with your local Sales/Service Representative. See the hardware manual for your instrument for information about installing these components or Table 1 for a comprehensive list of available components.


The following cover components can be replaced

- Right Side Panel-ISQ 7000
- Left Side Panel-ISQ
- Right Side Panel-ISQ Old
- Top Cover Panel-ISQ Systems
- ISQ Old Front Door Assembly (VPI)
- ISQ QD Front Door Assembly (No VPI)
- ISQ 7000 Front Door Assembly (No VPI)
- ISQ 7000 Front Door Assembly (with VPI)
- ISQ 7610 Front Door Assembly (no VPI)
- ISQ 7610 Front Door Assembly (with VPI)
- TSQ 8000 Evo Front Cover (with Vacuum Interlock Option)
- Right Side Panel-TSQ
- TSQ Duo Front Cover
- TSQ 9000 Front Cover (no VPI)
- TSQ 9000 Front Cover (with VPI)
- TSQ 9610 Front Cover (no VPI)
- TSQ 9610 Front Cover (with VPI)
- TSQ 9610 Front Cover (with NVAEI)
- Top Cover Panel-ISQ Systems
- Top Manifold Cover (glass)
- Front Door Hinge
- Front Door Hinge Support
- Front Door Latch
- Chassis Foot

Figure 12. Replaceable Covers



Right Side Panel-ISQ Old Quantity: Each System: ISQ, ISQ LT, ISQ QD System: ISQ, ISQ LT, ISQ QD ISQ QD • NOTE: To replace all the M4 x 10 mm screws on the right side panel, order four of PN 1R76913-0410. To replace the right side panel feet, order PN 1R3666-0207. Right Side Panel-ISQ 7000 Quantity: Each System: ISQ 7000



Thermo Scientific Part Number

1R120191-1020

NOTE: To replace all the M4 x 10 mm screws on the right side panel, order four of PN 1R76913-0410. To replace the right side panel feet, order PN 1R3666-0207.

1 Ordering Spare Parts Cover Components

Left Side Panel-ISQ Quantity: Each System: ISQ All System: ISQ All Image: Control of the system o

Thermo Scientific Part Number NOTE: The door you receive will be painted.



1R120407-0100

ISQ QD Front Door Assembly (No VPI) Quantity: Each System: ISQ QD System: ISQ QD Image: State of the system of the system: ISQ 7000 no VPI Isq 7000 Front Door Assembly (No VPI) Quantity: Each System: ISQ 7000 no VPI System: ISQ 7000 no VPI

Thermo Scientific Part Number

NOTE: The door you receive will be painted.





ISQ 7000 Front Door Assembly (with VPI)

Quantity: Each

System: ISQ 7000 with VPI

Thermo	Scientific	Part Number	

NOTE: The door you receive will be painted.

ISQ 7610 Front Door Assembly (with VPI)

Thermo Scientific Part Number

NOTE: The door you receive will be painted.

System: ISQ 7610 with VPI

1R120191-0110





Quantity: Each

1R120191-0100

ISQ 7610 Front Door Assembly (no VPI)	Quantity: Each
	System: ISQ 7610 without VPI
Thermo Scientific Part Number	1R120191-0210
NOTE: The door you receive will be painted.	
Top Cover Panel-ISQ Systems	Quantity: Each
	System: All ISQ systems
,	
Thermo Scientific Part Number	1R120412-0001

NOTE: To replace all the M4 x 10 mm screws on the top cover panel, order two of PN 1R76913-0410.



NOTE: To replace all the M4 x 10 mm screws on the right side panel, order four of PN 1R76913-0410. To replace the right side panel feet, order PN 1R3666-0207.

Thermo Scientific Part Number

1R120541-0200



Thermo Scientific Part Number

1R120192-0200

TSQ 9000 Front Cover (with VPI)	Quantity: Each
	System: TSQ 9000 with VPI
Thermo Scientific Part Number	1R120192-0100
TSQ 9610 Front Cover (no VPI)	Quantity: Each
	System: TSQ 9610 no VPI

Thermo Scientific Part Number

1R120192-0220

FSQ 9610 Front Cover (with VPI)	Quantity: Each
	System: TSQ 9610 with VPI
hermo Scientific Part Number	1R120192-0120
SQ 9610 Front Cover (with NVAEI)	Quantity: Each
	System: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-0121

Thermo Scientific

Left Hand Front Panel


Quantity: Each

*	•

Thermo Scientific Part Number

1R120541-0003

NOTE: To replace all the M4 x 10 mm screws on the right side panel, order four of PN 1R76913-0410.



Thermo Scientific Part Number

1R120541-0001

NOTE: To replace all the M4 x 10 mm screws on the top cover panel, order two of PN 1R76913-0410.

Top Manifold Cover (glass)	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120401-3000
Front Door Hinge	Quantity: Each
	System: All
Thermo Scientific Part Number	1R3434-5000

NOTE: To replace all the M3 x 8 mm (flathead) screws on a front door hinge, order two of PN 1R76937-0308. There are two hinges on the front door.



Front Door Hinge Support	Quantity: Each
Thermo Scientific Part Number	1R120445-0001

NOTE: To replace the M4 x 10 mm screw on the front door hinge supports, order PN 1R76913-0410.

Front Door Latch



Thermo Scientific Part Number

1R76483-3000

Quantity: Each

System: All

NOTE: To replace all the screws on the front door latch, order two of PN 1R120414-0010. To replace all the M3 locking nuts, order two of PN 1R76944-0100.

Chassis Foot	Quantity: Each
1	System: All
Thermo Scientific Part Number	1R3666-0206

NOTE: There are four feet on the bottom of the ISQ Series instruments, and six on the TSQ Series instruments, so to replace them, order four or six of PN 1R3666-0206. To replace the M4 x 12 mm screw on each chassis foot, order 1R76913-0412.

Pump Components

The following pump components are available for the ISQ and TSQ GC-MS instruments. Refer to the hardware manual for your instrument for installation instructions and to Table 1 for a complete list of pump components.

- Small Capacity Turbomolecular Pump (75 DX)
- Medium Turbomolecular Pump
- Large Turbomolecular Pump
- Convectron Gauge
- Rough Pump (RV3)
- Foreline Adapter

Small Capacity Turbomolecular Pump (75 DX)

Quantity: Each


System: All ISQ systems



Thermo Scientific Part Number

1R119268-0004

NOTE: If you are upgrading from the 70 DX pump, order the conversion kit PN 1R120567-0001.




Convectron Gauge	Quantity: Each
	System: All
Thermo Scientific Part Number	A0105-00501
Rough Pump (RV3)	Quantity: Each
	System: All
Thermo Scientific Part Number	1R76505-3007
Foreline Adapter	Quantity: Each
	System: All

MS System Tools

The MS systems ship with a toolkit that contains all the tools you will need to perform maintenance on the instrument. It does not contain any consumable items or cleaning supplies. You can purchase any item in the toolkit individually or you can order a complete kit. Contact your local Sales/Service Representative to place an order. See Table 1 for a comprehensive list of available components.

The following toolkit components can be replaced.

- MS Toolkit
- ISQ 7000, ISQ 7610, TSQ 9000, and TSQ 9610 VPI MS Toolkit
- NVAEI MS Toolkit
- VPI MS Toolkit
- Source Exchange Tool
- NeverVent AEI Source Exchange Tool
- Small Source Removal Tool
- NVAEI Small Source Removal Tool
- Bushing in the Source Exchange Tool
- Seal in the Source Exchange Tool
- Clip in the Source Exchange Tool
- Column Measuring Tool
- Source Holder
- T10 Torxhead Key
- T20 Torxhead Key
- T30 Torxhead Key
- Forceps, 8 in.
- Wrench, open-ended, 1/4-in., 5/16-in.
- Wrench, open-ended, 3/8 in., 7/16-in.
- Source Plug
- Source Plug Holder
- Source Plug O-ring



NVAEI MS Toolkit	Quantity: Each
	Systems: TSQ 9610 NVAEI Systems
Thermo Scientific Part Number	1R120189-6200
VPI MS Toolkit	Quantity: Each
	Systems: ISQ, ISQ LT, TSQ 8000, TSQ 8000 Evo


Thermo Scientific Part Number

1R120467-0002

Source Exchange Tool	Quantity: Each
	Systems: All VPI Systems
hermo Scientific Part Number	1R120406-2000
NeverVent AEI Source Exchange Tool	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-0003
Thermo Scientific Part Number Small Source Removal Tool	1R120189-0003 Quantity: Each

NVAEI Small Source Removal Tool	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-1017
Bushing in the Source Exchange Tool	Quantity: Each
	Systems: All VPI Systems
Thermo Scientific Part Number	1R120406-2203
Seal in the Source Exchange Tool	Quantity: Each
	Systems: All VPI Systems
Thermo Scientific Part Number	1R120406-2204

Clip in the Source Exchange Tool	Quantity: Pkg of 5
	Systems: All VPI Systems
Fhermo Scientific Part Number	1R76483-2102
Column Measuring Tool	Quantity: Each
	Systems: All
Thermo Scientific Part Number	1R120461-0010
	1R120461-0010 Quantity: Each
Thermo Scientific Part Number Source Holder	



Forceps, 8 in.	Quantity: Each
	Systems: All
Thermo Scientific Part Number	1R76360-0008
Wrench, open-ended, 1/4-in., 5/16-in.	Quantity: Each
3	Systems: All
Thermo Scientific Part Number	1R76360-0109
Wrench, open-ended, 3/8 in., 7/16-in.	Quantity: Each
3	Systems: All
Thermo Scientific Part Number	1R76360-0108
Source Plug	Quantity: Each
	Systems: ISQ 7000 and TSQ 9000 VPI systems


Source Plug Holder	Quantity: Each
	Systems: ISQ 7000 and TSQ 9000 VPI systems
Thermo Scientific Part Number	1R120589-1050
Source Plug O-ring	Quantity: Each
0	Systems: ISQ 7000 and TSQ 9000 VPI systems
Thermo Scientific Part Number	1R3816-202


Upgrade Equipment

The following upgrade equipment is available for the ISQ and TSQ GC-MS instruments. Contact your local Sales/Service Representative to place an order. Once you receive the equipment, refer to the hardware manual for your instrument for installation information or Table 1 for a comprehensive list of available components.

- CI Ion Source Cartridge
- AEI Ion Source Assembly
- CI Reagent Gas Flow Module
- CI Ion Volume
- EI/CI Ion Volume
- Hydrogen Ion Volume
- Dynode and Multiplier Power Supply and Cables (EI/CI)
- AEI Ion Volume
- AEI Ion Source Sleeve
- AEI Filament Retainer
- AEI Filament Spacer
- AEI Filament Shield
- AEI Ion Volume Insulator
- AEI Filament Insulator
- AEI Front Plate Assembly
- Source Magnet Retainer
- AEI Case Insert
- AEI Case
- AEI Magnet Retainer Bracket
- AEI Magnet
- AEI Source Magnet Assembly
- NeverVent AEI Ion Source Cartridge
- NeverVent AEI Ion Volume
- NeverVent AEI Filament
- NeverVent AEI Ion Cartridge Sleeve
- NeverVent AEI Filament Insulator
- NeverVent AEI Filament Cap
- NeverVent AEI Filament Spring

- Ion Gauge Option
- Ion Gauge Controller PCB
- Ion Gauge Tube Shield
- Large Turbomolecular Pump (nEXT 300D)
- Dust Filters
- Ion Gauge Mount
- Direct Insert Probe
- Direct Exposure Probe











Figure 15. NeverVent AEI Ion Source Cartridge Components R





AEI Ion Source Assembly	Quantity: Each
	Systems: ISQ 7000 and TSQ 9000
Thermo Scientific Part Number	1R120602-1300
CI Reagent Gas Flow Module	Quantity: Each
Thermo Scientific Part Number	1R23331-0092
CI Ion Volume	Quantity: Each
	Systems: All
Thermo Scientific Part Number	1R120404-4112

EI/CI Ion Volume	Quantity: Each
	Systems: All
Thermo Scientific Part Number	1R120404-4113
Hydrogen Ion Volume	Quantity: Each
	Systems: ISQ, ISQ QD, ISQ LT, TSQ 8000, TSQ 8000 Evo, TSQ Duo
Thermo Scientific Part Number	1R120404-4115
Dynode and Multiplier Power Supply and Cables (EI/CI)	Quantity: Each
	Systems: All

Thermo Scientific Part Number

1R120361-0007

NOTE: To replace all the M4 x 10 mm screws on the dynode and multiplier power supply, order two of PN 1R76913-0410. To replace the M4 x 16 mm screw, order PN 1R76913-0416.

AEI Ion Volume	Quantity: Each
	Systems: ISQ 7000, ISQ 7610, TSQ 900, TSQ 9610
Thermo Scientific Part Number	1R120602-0104
AEI Ion Source Sleeve	Quantity: Each
	Systems: ISQ 7000, ISQ 7610, TSQ 900, TSQ 9610
Thermo Scientific Part Number	1R120602-0101
AEI Filament Retainer	Quantity: Each
	Systems: ISQ 7000, ISQ 7610, TSQ 900, TSQ 9610
Thermo Scientific Part Number	1R120602-0102

AEI Filament Spacer	Quantity: Each
	Systems: ISQ 7000, ISQ 7610, TSQ 900, TSQ 9610
Thermo Scientific Part Number	1R120602-0103
AEI Filament Shield	Quantity: Each
•)	Systems: Systems: ISQ 7000, ISQ 7610, TSQ 900, TSQ 9610
Thermo Scientific Part Number	1R120602-0105
AEI Ion Volume Insulator	Quantity: Each
	Systems: ISQ 7000; TSQ 9000
Thermo Scientific Part Number	1R120602-0106
AEI Filament Insulator	Quantity: Each
3	Systems: ISQ 7000; TSQ 9000
Thermo Scientific Part Number	1R120602-0108


AEI Front Plate Assembly

Thermo Scientific Part Number

Quantity: Each

1R120191-0014

Systems: ISQ 7000; TSQ 9000 with AEI Source



AEI Case Insert	Quantity: Each
	Systems: ISQ 7000; TSQ 9000 with AEI source
Thermo Scientific Part Number	1R120602-1100
AEI Case	Quantity: Each
	Systems: ISQ 7000; TSQ 9000 with AEI source

Thermo Scientific Part Number

1R76524-0020

AEI Magnet Retainer Bracket	Quantity: Each
	Systems: ISQ 7000; TSQ 9000 with AEI source
Thermo Scientific Part Number	1R120602-1104
AEI Magnet	Quantity: Each
	Systems: ISQ 7000; TSQ 9000 with AEI source
Thermo Scientific Part Number	1R76253-2000
AEI Source Magnet Assembly	Quantity: Each
	Systems: ISQ 7000; TSQ 9000 with AEI source
Thermo Scientific Part Number	1R120602-0200
NeverVent AEI Ion Source Cartridge	Quantity: Each
	Systems: TSQ 9610 with NVAEI


NeverVent AEI Ion Volume	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-2001
NeverVent AEI Filament	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-6000
NeverVent AEI Ion Cartridge Sleeve	Quantity: Each
	Systems: TSQ 9610 with NVAEI

NeverVent AEI Filament Insulator	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-2036
NeverVent AEI Filament Cap	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-2035
NeverVent AEI Filament Spring	Quantity: Each
	Systems: TSQ 9610 with NVAEI
Thermo Scientific Part Number	1R120189-3002

Ion Gauge Option

Quantity: Each


System: All



Thermo Scientific Part Number

1R120560-0020


NOTE: To replace the o-ring on the ion gauge, order PN 1R3814-116. To replace the ion gauge tube, order PN A0105-06003.



on Gauge Tube Shield	Quantity: Each
	System: All
hermo Scientific Part Number	1R119605-0012
Large Turbomolecular Pump (nEXT 300D)	Quantity: Each
	System: All
hermo Scientific Part Number	1R119268-0002

Ion Gauge Mount


Quantity: Each



Thermo Scientific Part Number

1R120416-0002

Dust Filters	Quantity: Each
	System: All
Thermo Scientific Part Number	1R120442-1000



Note For information about ordering parts for the probes, see the *Direct Probe System User Guide*.

1 Ordering Spare Parts Upgrade Equipment